ABY-26008 CR.171 789

.

•

by

Pramode K. Bhagat
Principal Investigator

Prasad K. Kadaba Co-Investigator

Project Number - NAS9-15452 - Final Report

Project Period - November 1, 1977 - September 30, 1983

University of Kentucky
Wenner-Gren Biomedical Research Laboratory
Department of Mechanical Engineering
Lexington, Kentucky

# TABLE OF CONTENTS

| Acknowledgements  | i   |
|---|-----|
| Abstract  | v   |
| Introduction  | 1   |
| Background  | 1   |
| Plethysmography   | 1   |
| Relationship of Orthostatic Tolerance to Limb volume changes                  | 3   |
| Requirements for a Space Qualified Plethysmograph                             | 5   |
| Ultrasonic Plethysmography  | 6   |
| Single Chord Ultrasonic Plethysmograph  | 8   |
| Two Chord Ultrasonic Limb Volume System                                       | 31  |
| Ambulatory Ultrasonic Plethysmography   | 44  |
| Transducer Development  | 50  |
| Summary and Conclusions   | 66  |
| Theses, Dissertations and Publications  | 67  |
| Appendix A - Background Data on Transducers                                   | 70  |
| Appendix B - Test Procedure and Results Using Prototype ULP                   | 81  |
| Appendix C - Description of Individual Circuit Blocks of Two Transmitter ULVS | 108 |
| References  | 142 |

#### **ACKNOWLEDGEMENTS**

The development of ultrasonic plethysmography from a technical concept to a flight certified unit capable of measurements in space (zero G) environment has been the result of long term collaborative between the principal investigator and his associates with the National Aeronautic and Space Administration (NASA) personnel. This project was interdisciplinary in nature and required contributions from a variety of individuals.

NASA personnel with direct responsibility for technical monitoring of the project are gratefully acknowledged. In particular, Drs. Wayland Hull and Carter Alexander supported the initial development of the instrumentation detailed here. The continuous encouragement and long term support of Mr. Cletis Booher and his associates are sincerely appreciated.

The principal investigator is indebted to Dr. J. F. Lafferty, Director, Wenner-Gren Research Laboratory, for his continued support in terms of personnel, equipment and facilities for this project. The efforts of Mr. Tom Patrick, in the redesign and implementation of his pinger and receiver circuit electronics are gratefully acknowledged. It is a pleasure to acknowledge the assistance of Don Bowman, Charles Lewis, Dr. M. P. Kadaba, and Ms. Vicky Ruszcyzk in the development and evaluation of the single transmitter ULVS. Ms. Cindy Green assisted in the development of two transmitter transducer holder.

Also acknowledged is the support of Dr. W. T. Nickell, Mark Shafer, Dave Brewer and Ms. Julie Skidmore in the development and evaluation of two transmitter ULVS. The principal investigator is deeply indebted to the chief project engineer, Mr. V. C. Wu, for his participation in all phases of the reported work and his abilities to carry out all the needed tasks in a timely manner.

The assistance and support of the technical group at Denver Research Institute led by Whitney Patten towards the development of prototype flight certified ULVS is sincerely acknowledged.

### LIST OF FIGURES

- Figure 1 Schematic of Single Transmitter ULVS
- Figure 2 System Timing and D/A Convertor Schematic
- Figure 3 Avalanche Pulse Generator (Pinger) Circuit
- Figure 4 Ultrasonic Receiver Amplifier and Comparator
- Figure 5 Static Measurements of Dimensional Changes on the Lower Leg Using ULVS
- Figure 6 Cross Section View of Lower Limb at the Level of Plethysmograph Application
- Figure 7 Schema of System for Inducing Venous Occlusion in Lower Limbs
- Figure 8 Relative Placement of Ultrasonic Transducer and Whitney Strain Gauge
- Figure 9 Photograph Showing Actual Locations of Ultrasonic Transducers and Whitney Strain Gauge
- Figure 10 Calibration Curve for Whitney Strain Gage as Recorded on Chart
- Figure 11 Dynamic Comparison of Whitney Gauge and ULVS Measurements
- Figure 12 Computation of Arterial Flow and Compliance Level
- Figure 13 Correlation Between Initial Slope of Ultrasonic and Whitney Strain Gage Plethysmograph
- Figure 14 Average Initial Slopes (Ultrasound and Strain Gage) as a Function of Cuff Pressure
- Figure 15 Compliance Time Versus Cuff Pressure
- Figure 16 Computation of Cross-Sectional Area from Two Chord Length
  Measurements
- Figure 17 Two Transmitter ULVS

Recorder

- Figure 18 Block Diagram of ULVS Electronics
- Figure 19 System Electrode Placement
- Figure 20 Comparison of ULVS and Whitney Dimension Measurement Under Venous Occlusion
- Figure 21 Design of Flat and Hemispherical Transducers

- Figure 22 Construction of Receiver and Transmitter Elements
- Figure 23 The Developed Ultrasonic Transducers
- Figure 24 Off Axis Performance of Transducers
- Figure 25 Directional Characteristics of Ultrasonic Transducers
- Figure 26 Comparison of Performance of Different Transducer Combinations
- Figure 27 Domain Alignment in Ferroelectric Materials
- Figure 28 The Ultrasonic Field of Circular Disk Transducer
- Figure 29 Data Collection and Analysis Block Diagram
- Figure 30 Mechanical System
- Figure 31 ULP Transducer Output Pulse Shape
- Figure 32 Digital Output in EDI Mode (2400 Baud)
- Figure 33 Output Pulse with 50  $^{\Omega}$  Termination
- Figure 34 Amplitude Spectrum of the Output Pulse
- Figure 35 Pulse Shape with Transducer Load
- Figure 36 Tape Recorded Signal Played Back with 1K Load and Captured by Biomation A/D Converter
- Figure 37 Time Course of Dimension Changes Under 50 mm Hg Venous Occulsion Procedure
- Figure 38 Time Course of Limb Dimension Changes Under 500 Head Up Tilt
- Figure 39 ULP Data Played Back From the Recorded Venous Occlusion (50 mm Hg) Data
- Figure 40 Passive Monitoring of Data Through a Tilt Table Experiment
- Figure 41 Schematic of Processor and Memory for the Two Transmitter ULVS
- Figure 42 Schematic of Controller Logic, Transmitter Select and Counter for Two Transmitter ULVS
- Figure 43 Schematic of Power Supply Circuits Used in Two Transmitter ULVS
- Figure 44 Flow Graph of Software Implementation
- Figure 45 Subroutine "NOTYET"

# **ABSTRACT**

The mechanisms underlying the reductions in orthostatic tolerance with weightlessness are not well established. Contradictory results from measurements of leg volume changes suggest that altered venomotor tone and reduced blood flow may not be the only contributors to orthostatic intolerance. It is felt that accurate limb volume system which is insensitive to environmental factors will aid in better quantification of the hemodynamics of Of the various limb volume techniques presently available, ultrasonic limb volume system has proven to be the best choice. system as described herein is free from environmental effects, simple to operate and causes negligible radio frequency interference The segmental ultrasonic plethysmograph is expected problems. provide a better measurement of limb volume change since it is on cross sectional area measurements. Limb volume change is computed from the two cross-sectional areas measured at different sites length apart) using an equivalent right circular cylinder geometry. A space qualified version of our device has been built by Research Institute and is expected to be used in a future shuttle mission.

#### I. INTRODUCTION

The study and evaluation of cardiovascular deconditioning during space flight impose unique requirements on the limb volume measuring system (LVMS). The LVMS must operate, in both earth and zero gravity environments, be simple to calibrate and respond to both small and large changes in limb volume and be largely insensitive to humidity, ambient temperature and pressure variations.

The objective of this project was to develop a new technique for real time measurement of lower limb volume changes in space. The technique developed in this study, the ultrasonic plethysmograph, is based on transit time measurement of an acoustic pulse propagating through tissue. A microprocessor based ultrasonic lower limb volume measuring system was developed as a prototype system for use in space flight applications. It is based on the assumption that the cross section of the lower limb can be approximated by a circle. On this basis two independent chord length measurements using ultrasonic through transmission technique are made at a given site on the lower limb. This data is used by the microprocessor to compute the area of a circle passing through the measured chord lengths.

This report details the development of the ultrasonic limb volume measuring system in to a flight certified unit which is expected to be used in a future space flight mission.

### II. Background

### A. Plethysmography

Plethysmography is the recording of variations in the volume of a liquid or undissolved gas in biological tissue. Studies of

peripheral circulation have been carried out through a variety of methods which record/measure the volume of blood content in a limb (plethysmography).

principle of plethysmography was first introduced The by limb Glissonio (1677) who used water displacement for measurement (1). Variations of this method, some using air as the displaced fluid, are still used for the measurement of arterial flow into segments or distal portions of limbs. Another frequently used plethysmograph is the Whitney strain gage, introduced in 1953 (2). The Whitney, or mercury-in-rubber, strain gage is attached around the segment in question and any relative change in the electrical resistance of the strain gage is used as a measure of a change circumference of that limb segment. The accuracy of this method plethysmography was compared to that of water displacement plethysmography by Paulev, et al.(1). They found that "the accuracy of the water and strain-gage plethysmograph is identical, and the two instruments closely follow random variations in blood flow." However, they further state that "The results of the strain-gage plethysmograph are dependent on the actual location of the strain gage."

plethysmography include impedance Other methods of capacitance plethysmography. Impedance plethysmography (rheography) makes use of the fact that the electrical impedance of a limb segment varies with its blood content. Capacitance plethysomgraphy involves the use of the skin as one plate of a capacitor. An electrode to a band around the limb segment completes the capacitor. attached Assuming a purely radial expansion with increased limb volume, in capacitance between the skin and the electrode change is proportional to the change in volume of the limb segment. These and

other currently used methods of plethysmography are discussed more thoroughly by Sigdell (3). A comparative analysis of their feasibility for space flight applications may be found in an earlier report from this laboratory (4).

B. Relationship of Orthostatic Tolerance to Limb Volume Changes

Reductions in orthostatic tolerance due to exposure to zero G environment were first reported during the late flights of Merucry when tilt table tests conduced post-flight revealed moderate orthostatic hypotension. These findings were also observed during the Gemini program (flight duration 3-14 days). However, responses to the 70 degree upright tilt usually returned to normal within fifty hours after splash down (5,6).

During the Apollo program (7), a new test procedure, lower body negative pressure (LBNP) was used to study orthostatic tolerance. LBNP has several advantages over tilt table testing, some of which are:

- 1. No movement of subject is required which leads to simplified instrumentation and stable physiological signals.
- 2. Ease of application of controlled stress (magnitude and duration) to the cardiovascular system.
- Suitability of use in the weightless condition.

LBNP studies were performed for most Apollo crews, and the response of cardiovascular system during weightlessness was inferred from pre-flight and post-flight data. The results demonstrated reduced orthostatic tolerance, which returned to normal within 48 hours. Supine measurements of maximal calf circumference were also performed, before and after flight, of 24 crewmen. During post-flight evaluation 67% of crewmembers showed significantly reduced calf circumference, and this was not totally regained approximately 120

hours after splash down (7,8).

During the Skylab Program, the LBNP experiments were performed during missions 2,3 and 4 to determine the extent and time course changes in orthostatic tolerance during weightlessness (9). An additional objective was to determine whether the inflight data can be used for prediction of post-flight orthostatic tolerance status. LBNP produced exaggerated blood pressure and heart rate responses during the first inflight test of Skylab 2 crew, but did not show a trend towards the pre-flight levels during the 28-day flight. observed that cardiovascular responses to LBNP tended to become more stable after six to eight weeks of flights, but no trend toward preflight levels could be seen. Johnson et al (9) suggest that total blood volume may already have been reduced by the time first inflight LBNP tests were performed, and that a significant fraction of fluids previously located in the legs was now accommodated by veins intersitial tissues of the upper part of the body. Thornton Hoffler (10) conducted studies of lower limb during the Skylab 4 mission to test the hypotheses that changes in hemodynamics of legs increased blood pooling and reduction in cardiac output is the with most probable cause of orthostatic intolerance following spaceflight. In spite of the individual variability, the data showed increased blood flow above pre-flight and post-flight levels. compliance changes were high for both the scientist pilot and pilot whereas for the commander there was little change, and post-flight data showed return to pre-flight values of compliance. Based on scant data from three crew members, the authors were unable to further prove their hypothesis. However, their findings support the need for

careful assessment of pressure flow relationship in the lower limb.

Measurement techniques used in these redistribution studies include estimates of leg volue using multiple tape measurements of limb circumference, standard and infrared photography, and center of mass measurement (7,8). During the Skylab 3 and 4 missions, volume measurements were made using a space qualified version of capacitance limb volume system developed by Barnett et al (11). This system represented the best available at that time for measurement of limb volume, although it is affected by environmental factors, such as humidity and temperature variations. In addition, all of these techniques suffer from the difficulty that considerable crew effort is required to obtain data, or that the transducers are bulky and cannot be used during launch or normal activities. For this reason, no experiments could be performed during the first hours and days of zero-gravity.

- C. Requirements for a Space Flight Qualified Plethysmograph Based on previous research (4) it is felt that an ideal space flight plethysmograph should have the following characteristics:
  - 1. Operation in both zero and earth gravity environment.
  - 2. Insensitive to humidity.
  - 3. Insensitive to ambient temperature and pressure variation.
  - 4. Compatible in size and operation with LBNP test procedures.
  - 5. Simplicity of operation and calibration.
  - 6. Reliable and safe.
  - 7. Nonconfining and non-restrictive of limb dimensional changes.
  - 8. Insensitive to limb movement.
  - 9. Must respond to both small and large changes in limb volume.

### III. ULTRASONIC PLETHYSMOGRAPHY

Ultrasonic plethysmography is based on measurement of transit time required for an ultrasonic pulse to propagate from one given point to another. If the velocity of the propagating sound wave in the medium under study is known, then measurement of transit time leads to computation of the pulse traversed distance, d:

$$d = c.t \tag{1}$$

where c = velocity of sound in the medium (m/sec)

t = measured transit time (sec)

Velocity of propagation of acoustic waves have been measured and reported in selected animal tissues using pulse echo techniques. Goss et al (12) provide a comprehensive listing of published results on acoustic properties of mammalian tissues. For a given frequency and temperature, the velocity of propagation appears to be constant and the range of variation is quite small for soft tissues. The average velocity of sound in soft tissues (excluding fat tissue) is  $1560 \pm 20$  m/sec at 2.25 MHz (Bhagat et al (12), Kadaba (13)). The temperature dependence of sound velocity in the range of 20-40 degrees C shows a variation of less than 2m/sec/degree C. Since the temperature within the living tissue is maintained constant by the body within very narrow limits, temperature related errors can be safely ignored.

The use of ultrasonic transit time measurements for biological studies was first reported by Rushmer et al (14) who monitored left ventricular dimensions in open chested laboratory animals. Stegall and colleagues (15,16) have reported the instrumentation and results associated with implementation of the transit time principles.

Stegall et al (17) in an abstract dealing with noninvasive diagnosis of peripheral vascular disease suggested the use of an

ultrasonic dimension gauge to measure changes in limb lower However, we have been unable to find any details on its implementation. We were the first to demonstrate the actual feasibility of ultrasonic dimension gauge for noninvasive monitoring of lower limb dimensions (18). In this approach, two ultrasonic transducers (resonant at 2.25 MHz) are placed across the limb segment under study. One of the transducers is shock excited which causes it to vibrate at its natural frequency. This ultrasonic vibration is transmitted through the limb and is converted into an electrical signal by the receiving transducer. The time period between initiation of the shock pulse and its reception by the receiving transducer can, therefore, be accurately measured. With the knowledge of acoustic velocity in the lower limb, the segment dimension is easily obtained.

The most critical choice for ultrasonic plethysmography of lower limb segment is that of transducers' natural frequency. Ultrasonic waves are attenuated when propagating through biologic media. The amplitude of the propagating wave decreases by an exponential factor given by

$$P_{o}(f) = e^{-2\alpha(f)x} P_{i} (f)$$
 (2)

where  $P_{O}(f)$  is the amplitude at a distance x in the medium

 $\alpha(f)$  is the attenuation coefficient of the medium (dB/cm)

P; (f) is the input wave

and f is the frequency of the wave.

Measurements in our laboratory (19) indicate that the attenuation coefficient  $\alpha$  (f) is approximately linearly related to ultrasonic frequency ( $\alpha$  (f) = af<sup>b</sup>; b = 1.06) for excised muscle tissue in

experimental animals. In light of this it is imperative that the ultrasonic frequency be kept as low as possible.

As the ultrasonic frequency is lowered two constraints come into play. In most transit time measurement systems arrival of the ultrasonic pulse is threshold detected with an amplitude comparator. Thus the maximum error due to threshold adjustment can be  $\lambda/4$ , where  $\lambda$  is the wavelength of ultrasound. Assuming an average velocity of 1,560 m/sec, and calf length of 15 cms, in the range of 2-3 MHz the maximum threshold error is approximately 0.2 mm.

Further as ultrasonic frequency is lowered the thickness of crystal increases. Higher shock voltages are required at lower frequencies to set the transducer into vibration, due to its increased inertia. While this is not critical in most applications since the excitation voltages are usually in the range of 100-200v and are of short duration (4 ns), the authors feel that from a safety viewpoint every attempt should be made to keep these voltages lower.

We, therefore, chose frequency range of 2-3 MHz for the plethysmographic applications.

### SINGLE CHORD ULTRASONIC PLETHYSMOGRAPH

The single chord ultrasonic plethysmograph is schematically shown in Figure 1. The major components of this instrument are a pulse generator (Pinger) and receiver amplifier, a voltage comparator, a sync pulse generator and an elapsed time counter. The system is housed in a single unit, with the exception of the transmitting and receiving piezo-electric crystals. The crystals are cut from 1" square LZT-5 stock, resonant at 2 MHz, with approximately 3 mm diameter, and then wires are attached. The crystals are dipped in




FIGURE 1
SCHEMATIC OF SINGLE TRANSMITTER ULVS

clear epoxy to maintain the integrity of the transducer material over a period of time. The crystals are positioned on a limb such that the ultrasonic pulse will travel from the transmit crystal to the receiver crystal. The wave propagation time between crystals is measured by counting the number of cycles generated by a time base oscillator. Output from the system is available as a 12-bit binary number or an analog signal derived from the binary number through a digital-to-analog converter (DAC).

The measurement cycle begins on the leading edge of a 40 ns wide negative-going pulse. This pulse is generated by a monostable multivibrator (IC6, Fig. 2) triggered by a crystal controlled oscillator and frequency divider at a repetition rate of 1 KHz. This pulse induces a 150-V signal across a piezoelectric crystal. Due to this shock excitation, the crystal vibrates at its natural frequency, 2 MHz, resulting in the transmission of an ultrasonic signal across the limb.

The pinger circuit (Figure 3) works as follows:

When the trigger level is high (+15V) the base to emitter junction in transistor  $Q_1$ , is non conductive, therefore, no current flows through it. This also implies that transistor  $Q_2$  will also be off and no signal is available across the transformer to the ultrasonic crystal. However, when the trigger level is lowered to zero,  $Q_1$  and therefore  $Q_2$  are turned on. With  $Q_2$  in the operating mode a 200V signal is induced through the 1:1 transformer to excite the transmitting transducer. The transformer is included in the pinger circuit to provide isolation of 200V source from the test media (usually human subjects). The 200V signal is derived from a commercially available module (ERG) which is a dc to dc invertor with





FIGURE 2. SYSTEM TEMING AND D/A CONVEXTOR SCHEMMIC



AVALANCHE PULSE GENERATOR (PINGER) CIRCUIT


FIGURE 3

input of 15V.

Simultaneous to the triggering of the transmitter circuit, the elapsed time counter (IC's 10-13, Fig. 2) is cleared and the time base signal (32 MHz) is gated to the input (IC 10) of the counter. The count accumulates until the voltage comparator senses a signal from the ultrasonic receiver which exceeds the operator controlled threshold voltage.

The input signal to the ultrasonic receiver (Fig. 4) comes from a piezoelectric receiving crystal attached to the limb. the receiving crystal is connected to the primary winding to  $T_2$  which is balanced by a 1 kohm potentiometer and a 220 µf capacitor to provide common mode rejection. Transformer  $T_2$  also provides electrical isolation for the subject. The signal induced in the secondary of  $T_2$  is then amplified by a factor of 1,000 with two cascaded RF amplifier (LM371). The resulting signal serves as input to the voltage comparator (LM306).

Once the comparator senses the received signal, the time base signal being input to the counter is halted. After a short delay, allowing for propagation delays in the counter, the final count is latched in a separate register. The digital-to-analog converter changes the 12-bit register output to a proportional voltage which is recorded by an analog recorder. The register output can also be input into a computer system for data storage and analysis. A portion of the output can be suppressed by dialing in a number on a set of thumbwheel switches. The negative value of this number is loaded into the counter at the beginning of each measurement cycle thereby subtracting the number from every measurement. This feature allows



the operator to record only the variations from a baseline reading. The distance resolution of this system using a 32 MHz clock is given by resolution = (sound velocity/clock frequency) = 0.049 mm. This system has been described in literature (18).


Extensive experimentation was carried out using human subjects to investigate the reliability, accuracy and repeatability of the developed instrument. Experimental protocols are defined below:

## a) Static tests

With the subject in a supine position on a tilt-table, the instrument was tested by inducing a small precise changes in leg A 2 MHz transducer was mounted on each tip of a 5 diameter. With the micrometer rigidly clamped, the leg of the subject was positioned so the micrometer spanned the major chord length from near the crest of the tibia to the back of the qastroenemius. The micrometer head was advanced until the crystals made firm contact with the skin. The chord length was then varied by advancing the micrometer head in steps of 0.025 in (0.635 The change in the chord length as measured by the ultrasonic system is shown in Figure 5. As can be seen from the figure the system response discrete step changes from 11.74 cm diameter to 11.30 cm diameter is repeatable as well as accurate.

### b) Dynamic tests

A comparative study using Whitney Strain Gauge, (WSG), "Mercury in rubber" plethysmograph, was carried out to establish the dynamic performance of the system. If the limb segment involved in the plethysmographic measurement is assumed to be circular and the path between the ultrasonic transducers is assumed to be its diameter, then the cross-section area or volume per unit length of the limb




CHORD LENGTH ---

segment is directly proportional to the measured distance. Any change in volume is, therefore, reflected in a change in travel time of the ultrasonic pulse. The Whitney Strain Gauge measures changes limb. circumference of the Based on this hypothesis plethysmographic measurements by the WSG and the ultrasonic plethysmograph, placed over the same limb segment, should be linearly related.

A typical cross-section of the lower leg in a healthy adult, at the level of limb segment normally involved plethysmography studies, is shown in Fig. 6 with approximate transducer locations. As can be seen from the figure, the ultrasonic distance measured is more like a chord of the limb segment. However, the ultrasonic output and the Whitney Strain Gauge output should still be related by a constant.


In this experiment, venous occlusion was used to produce a discrete change in limb volume. A pressure cuff was applied to the proximal portion of the limb and the pressure in the cuff increased above the venous pressure but below the diastolic arterial pressure (usually between 20-60 mm Hg). Because the arterial of blood was not stopped, there would be a net increase in the volume of the limb below the cuff due to a pooling of blood in the segment. The rate of volume change with respect to time immediately after occlusion can be used to determine the rate of blood flow in the The volume will stabilize when the pressure in the veins arteries. reaches the cuff pressure. This is known as the compliance level and, along with the initial rate of volume change, is one of the important parameters in venous occlusion plethysmography. A more thorough discussion of the mechanisms of venous occlusion plethysmography is



SECTION 96

- (->) Approximate position of transducers
- 1. Tibia
- 2. M. flexor digitorum longus:
- 3. V. saphena magna et n. saphenus
- 4. V. et a. tibialis posterio
- 5. N. tibialis
- 6. M. soleus
- 7. Tendo m. plantaris
- 8. M. gastrocnemius
- 9. N. cutaneus surae medialis
- 10. V. saphena parva
- 11. N. cutaneus surae lateralis
- 12. M. gastrocnemius
- 13. A. et v. peronaea
- 14. Fibula
- 15. N peronaeus superficialis
- 16. Mm. peronaei longus et brevis
- 17. N. peronaeus profundus a. et v. tibialis anterior
- 18. Mm. extensores longi digitorum et hallucis
- 19. Membrana interossea cruris
- 20. M. tibialis anterior

FIGURE 6 Cross sectional view of lower limb at the level of plethysmograph application.



given by Sigdell (3).

The apparatus used for this study includes two ultrasonic transducers, a Whitney strain gauge, a system for obtaining venous occlusion of the lower limbs, and the single chord ultrasonic plethysmograph.

The venous occlusion system, shown in Figure 7, consists of a twelve gallon fixed displacement air compressor (a), a pressure attenuator (b) for limiting the pressure in a small pressurized tank (c), two normally closed solenoids (d), and standard thigh cuffs equipped with sphygmomanometers. Upon activation of the first solenoid in the system, the pressure in the thigh cuffs is increased in approximately a stepwise matter. The amount of pressure in each cuff is determined from the sphygmomanometers. Once the first solenoid is deactivated the pressure in the thigh cuffs remains constant until the second, release, solenoid is activated.

With the subject relaxed in a supine position, the ultrasonic transducers are positioned on the lower leg at the level shown in Figure 6, using the oscilloscope to determine the position corresponding to maximum signal. Once the optimum position is determined, the transducers are taped in place using the smallest piece of tape necessary for secure placement. The Whitney strain gauge is then placed over the same limb segment as shown in Fig. 8 schematically and a photograph is shown in Fig. 9. The leg is, then, elevated slightly to eliminate pressure on the calf muscle and to prevent drag on the strain gauge. The thigh cuffs are applied to both legs above the knee.

After a stable base line is obtained on both channels of the chart recorder, the initial readings of the ultrasonic digital display

) Normally closed ASCO solenoids electrically controlled by a two-way switch

12 gallon fixed displacement

air compressor

b) Pressure attenuator


Small pressure tank

૽

e) Standard thigh cuffs with sphygmomanometers

FIGURE 7

Schema of system for inducing venous occlusion in lower limbs.



- a. Ultrasonic transducer taped in place
- b. Whitney strain gage
- c. Ultrasound connection wires -- to oscilloscope, digital display, and chart recorder
- d. Strain gage connection wires--to gage display and chart recorder

# FIGURE 38

Relative placement of ultrasonic transducer and Whitney strain gage




FIGURE 9

and the strain gauge display are recorded. The thigh cuff pressure is then increased to 20 mm Hg and held constant until the chart recorder graph reaches a maximum (compliance level). Readings are again recorded for the strain gauge display and the ultrasonic digital display. The pressure is then released and the system allowed to stabilize. This procedure is repeated for pressures of 30, 40 and 50 mm Hg. In this project at least three complete sets of data (i.e. pressures of 20-50 mm Hg) were taken during each session. Data sets were taken for each leg under the conditions of bilateral and unilateral venous occlusion. However, the equipment did not permit collection of data from both legs simultaneously.

To facilitate the interpretation of chart recorder data in terms of physical parameters (i.e. change in length per unit time and percent change in length at compliance), it was necessary to perform calibration experiments for the Whitney strain gauge and the ultrasonic system. The calibration of the ultrasonic system has been described under static test.

The calibration procedure for the Whitney strain gauge was as follows: One end of the strain gauge was fixed and the other end attached to a vernier device. The system was adjusted to an initial position with the strain gauge straight but not extended. The strain on the strain gauge was then incremented in 1 mm steps to a total of 10 mm strain. The gauge reading was recorded at each increment to check the linearity of the readings. The strain gauge was then returned to its initial length in steps of 1 mm and the readings were again recorded. Figure 10 shows the static data for Whitney strain gauge.

A typical venous occlusion record is shown in Figure 11. Cuff




FIGURE 10 Calibration curve for Whitney strain gage as recorded on chart recorder




FIGURE 11: DYNAMIC COMPARISON OF WHITNEY GAUGE AND ULVS MEASUREMENTS

occluding pressure of 50 mm Hg was applied to the subject identified by pressure on arrow on the graph. The occlusion was maintained for a time period of approximately 2 minutes when the ultrasonic plethysmograph output reached steady state. The occluding pressure was released then and the record indicates a rapid emptying of the pooled blood followed by a slower decline to baseline values. Calculations of maximal blood flow and compliance levels can be made and the procedure is illustrated in Figure 12.


Figure 13 shows initial slope of the ultrasonic plethysmograph plotted against initial slope obtained using the Whitney strain gauge. The equation of the line shown in this figure, as determined by the least square error method, is

$$W = 2.235 u + 0.176$$

The above equation was derived using data from 30 independent experiment runs. The correlation is quite good (r=0.9077).

Figure 14 shows the average initial slope as a function of cuff pressure. As can be seen, the initial slope remains approximately constant for pressures of 20-40 mm Hg and drops off slightly when the cuff pressure is 50 mm Hg. One possible explanation is that 50 mm Hg occluding pressure may have caused occlusion of some arterial vessels and thus decreasing the inflow rate.

The compliance level showed a direct dependence on the applied cuff pressure. Both the ultrasonic plethsmograph and the Whitney strain gauge outputs leveled off around the same time frame. Figure 15 shows a plot of compliance time versus cuff pressure for two subjects. As expected time to reach final compliance level increased with the applied occluding pressures.



Time \_\_\_\_\_

Arterial inflow =  $\frac{dv}{dt}$ Compliance level =  $V_{C}$ 

Figure 12 Computation of arterial flow and compliance level




FIGURE 13: Correlation between initial slope of ultrasonic and Whitney strain gage plethysmographs




FIGURE 14: Average initial slopes (ultrasound and strain gage) as a function of cuff pressure




FIGURE 15: Compliance time versus cuff pressure

#### TWO CHORD ULTRASONIC LIMB VOLUME SYSTEM

single dimension measurement is not a true indicator of volume change in the lower limb. It was felt that a more accurate measure of limb volume could be obtained through computation of a cross-sectional area from two independent chord length measurements. independent chord lengths are measured ultrasonically The placement of ultrasonic transducers around the circumference of the We exploit the well known concept from analytic geometry that limb. points in a plane can define a unique circle. three implementation of this concept uses three transducers, located on the Two of these transducers, kept a fixed distance apart and limb. one side of the limb act as transmitters of ultrasonic energy. This energy is received by the third transducer which is located on the opposite side of the limb. Thus, two independent transit time measurements can be made at any given site on the With the knowledge of ultrasonic propagation velocity, the chord lengths can be computed as described earlier in this report. It may be noted that this computation also provides a sort of "averaging" of the ultrasonic chord length measurements, whenever the tissue expansion is approximately isotropic even when the cross-section noncircular. A microprocessor was incorporated into the system to provide real time computations of limb cross-section, as well as control the entire measurement process. Further the use of microprocessor makes it easy to change the computational and control routines thus facilitating subsequent refinements of the design.

#### Theoretical details:

The following mathematical scheme defines the procedure for computation of the cross-sectional area from chord length

#### measurements:

Consider Figure 16 where OA and OB are the chord length obtained from the measured transit times using transmitters (located at A and B respectively) and receiver transducer (located at 0). The distance, AB, between the transmitters is fixed and known =  $D_3$ cm and

$$OA = t_{OA}$$

$$OB = t_{OB}$$

Where  $t_{OA}$  and  $t_{OB}$  are times required for the ultrasonic pulse to be received at 0 from the sites at A and B respectively.

Without any loss of generality we can assumed that OA defines the x axis. Hence the x-coordinate of the center

$$x_C = D_1/2 \tag{1}$$

The equation of the circle passing through O, A and B with center  $(x_C, y_C)$  is given by:

$$(x-x_c)^2 + (y-y_c)^2 = r^2$$
 (2)

Where 'r' is the radius of the circle.

At the point O(0,0) the equation is

$$(D_1/2)^2 + y_c^2 = r^2$$
 (3)

and at the point B  $(x_2, y_2)$ 

$$(x_2-D_1/2)^2 + (y_2-y_c)^2 = r^2$$
 (4)

also

$$D_3^2 = D_1^2 + D_2^2 - 2D_1D_2 \cos\theta \tag{5}$$

solving equations 3, 4 and 5 for radius, r, we obtain




FIGURE 16. Computation of cross-sectional area from two chord length measurements

$$r^{2} = x_{c}^{2} + y_{c}^{2}$$

$$= D_{1}^{2}D_{2}^{2}D_{3}^{2}/(D_{1}+D_{2}-D_{3}) (D_{1}+D_{2}+D_{3}) (D_{1}+D_{3}-D_{2}) (D_{1}-D_{2}-D_{3})$$
(6)

Area of the circle passing through O, A and B is  $_{\pi}$   $_{r}^{2}$ 

which can be easily computed with knowledge or radius, r.

Equation 6 was simulated on both a minicomputer and a microcomputer for computational accuracy using several chord length sizes. All calculations were carried out in floating point arithmetic. program, written in assembly language using Intel 8080 microprocessor chip, took approximately 1.5 k memory space. The accuracy of microprocessor computations was compared with that obtained using a DEC PDP-11/10 minicomputer and the results were in excellent agreement with each other. The truncation error in area computations using  $D_1$  = 6.54 cm;  $D_2 = 7.55$  cm and  $D_3 = 1.85$  cm was  $5*10^{-3}$ %, an insignificant should be noted that these computations are affected by error. accuracy in measurement of D3 (which must be accurate to within 0.1 large D<sub>3</sub> provides very accurate computations but leads practical problems in implementation of the concept since the bony tissue, tibia, must not be in path of ultrasonic transmission. addition, operational difficulties with actual placement of transducers dictate the transmitters be mounted on the same assembly. We have determined that D3 of the order of 2.5 cm provides a good compromise between computational accuracy requirements and operational ease. (21).

Details of the developed instrument:

# A. <u>Hardware Description</u>:

The developed instrument (Fig. 17) is housed in a Tektronix TM An SC 502 oscilloscope is included to monitor quality of received signals. The instrument is modular transmitters, receiver, and processor circuits housed in separate Transmitter and receiver circuits are identical described earlier for chord the single ultrasonic plethysmograph. Two separate transmitter (pinger) circuits Either of these pinger circuits can be selected by processor, providing independent measurement of the two chord lengths. The piezoelectric transducer elements, which are mounted in speciallymolded housings for attachment to the limb, are connected to the instrument by unshielded twisted pair cables.

Processor and Memory: A block diagram of the implemented circuits is shown on Figure 18. The system uses one Intel 8080 CPU with a clock generator (8224), which provides 18 MHz for the counter and the timing frequency of 2 MHz for the CPU. 256 byte random access memory (RAM) is used for storage of generated data and a 2kb read only memory (ROM) is available for storage of software codes and constants needed for computation of cross-sectional area. The data lines from CPU are latched through a data bus and control status bus driver in order to increase the drive capacity of the CPU. Two peripheral inteface adapters (PIA) are used: one connects the data bus to LED display through four BCD to seven segment decoder/drivers while the other interfaces the counter to the CPU. Due to the available number of lines on the PIA chosen the output display is limited to less than 200 cm<sup>2</sup> (The 5th digit of the display is wired so that it will show





FIGURE 17: Two transmitter ULVS



PIGURE 18. BLOCK DIAGRAM OF ULUS ELECTRONICS.

blank or "1" only). This restricts the maximum calf diameter to less than 16.5 cms. For the normal adult population this is not a serious limitation. The power for all the units is derived from the TM 515 power module with suitable regulators to obtain  $\pm$  12 V and  $\pm$  5 V.

Transmitter Select and Counter Control Logic: The 12-bit binary counter is interfaced with the CPU through PIA as shown. The 18-MHz clock frequency provides 56 ns resolution for time measurements. Following transmitter activation, the comparator output is inhibited for 10  $\mu$ s to allow electrical interference from the transmitting pulse to decay. This block also contains transmitter select and level conversion circuitry necessary for producing a pulse from one of the two transmitters. The command pulse which fires the selected transmitter also starts the counter. The counter continues to count until it is stopped by a level change signal from comparator.

A three-position switch (mode switch) is used to select the modes of operation of the unit. Each of the first two positions causes one of the two transmitters to be "pinged" at a 2-KHz rate. The receiver amplifier and comparator outputs are displayed on the oscilloscope so that the quality of the received signal can be assessed. The pinger and receiver amplifier circuits are shown in Figures 3 and 4 of this report.

Transducer Assembly: The two transmitting crystals are mounted in a plexiglass assembly which is curved to fit the shape of the leg. The distance between the two crystals is set at approximately 2.5 cm. This spacing was selected as a compromise between computational accuracy and ease of fabrication. Crystals are cut from sheets of lead zirconium titanate transducer material (LTZ-2), with a resonant

frequency of approximately 2 MHz. The transmitting crystals are relatively large (5 mm diameter), in order to ensure an adequate signal level at the receiver. The receiving crystal is smaller (2 mm). Use of a small crystal makes the receiver less directional, but does not seriously reduce signal amplitude. All the crystals are fitted with lenses which produce a wider beam width than would be provided by the crystal alone.

# B. <u>Software Description</u>:

The instructions necessary to operate the instrument are contained within the 2K ROM. This ROM contains, in addition to operating instructions, subroutines for floating point arithmetic and for conversion of floating point numbers into a code suitable for output to the LED display. Computations are carried to five significant figures.

The operating routines are designed to make use of the instrument as simple as possible. When the three-position selector switch is in the OPERATE position, pushing RESET causes the instrument to output a series of voltage levels on the analog output. These voltage levels produce a calibration record for the strip chart recorder.

The instrument next measures two chords and computes a baseline value for the cross-sectional area. Chords are then measured at a rate of about 10/s. Simultaneous with the analog output, the computed cross-sectional area is sent to the display. The analog output is scaled so that full scale deflection of the recorder represents +6% to -2% volume change.

Two digital filters are used to improve the appearance of the analog output. A slew rate filter serves to limit the effect of

occasional false comparator triggers. This filter functions by limiting to three steps the maximum change allowed in the 8-bit D/A convertor in one sampling interval. Thus, very large and sudden deviations are not passed to the output.

A one pole low-pass filter is used to reduce the effect of digitization and other noise on the analog output. This filter requires storage of one previous value of the output and one multiplication. Implementation, therefore, requires only a minimum of processor time and memory. In the current version of the software, the time constant of the filter is set at about 3 s.

Additional details regarding actual hardware and software implementation can be founded in Appendix C.


# C. Operation of the Instrument:

To use the instrument, ultrasonic transducers must be properly attached to the subject's leg. Transducers are mounted using exercise EKG adhesive collars (Beckman). The transducers are attached first at midcalf level in such a way as to avoid passing the ultrasonic pulse through the bones of the leg. This is conveniently accomplished by placing the received transducer just to the inside of the tibia. The transmitter assembly is then moved around the rear half of the calf until adequate signals are obtained. By use of the three-position selector switch, the signal from either of the two transmitting crystals can be displayed on the oscilloscope screen. The position of the receiving transducer is adjusted until an acceptable signal is obtained from both transmitters. When such a position is found, the selector switch is set to the OPERATE position. Pushing RESET causes instrument to produce a series of calibration levels before beginning to measure and output fractional changes in leg volume.

analog output of the ULVS should be connected to a strip-chart recorder with  $\pm 2.5$  V full scale sensitivity.

Comparison of ULVS with other plethysmographic techniques:


The developed two transmitter ULVS was used in a study by NASA (22) to compare the sensitivity and accuracy of our methodology. A population of thirty subjects (50% male) was used. These subjects were non-paid volunteers (age 18-45), in good health with no evidence of varicosity. Tests were carried out in an environmental chamber with controlled temperature and humidity. The Whitney Strain Gauge was used as a standard and compared against our ULVS and an impedance plethysmograph. Figure 19 shows the system electrode placements and the results of a 50 mm Hg venous occlusion are shown in Figures 20. The correlation coefficient between the WSG and ULVS for the computed limb volume was 0.996. The maximum volume indicated by the ULVS was consistently higher than the corresponding WGS output. The impedance plethysmograph gave similar results when compared against WSG (22).



SITE OF ONE ULTRASOUND ELECTRODE IN WHITNEY/ULTRASOUND TESTING, SECOND ELECTRODE APPROXIMATELY 135° AROUND BELLY OF GASTROCNEMIUS MUSCLE.

= SITE OF IMPEDANCE ELECTRODES IN WHITNEY/IMPEDANCE TESTING.

TOTAL TO CYCTEM ELECTRONE DI ACEMENT



COMPARISON OF ULVS AND WHITNEY DIMENSION MEASUREMENT UNDER VENOUS OCCLUSION FIGURE 20.

#### AMBULATORY ULTRASONIC PLETHYSMOGRAPHY:

The following technical requirements were developed in consultation with responsible NASA personnel for fabrication of an ambulatory ultrasonic limb plethysmotraph, ULP.

# Overall System Requirement:

The contractor (equipment fabricator) shall be required to produce a flight qualified, portable four transmitter ultrasonic measuring sections plethysmograph capable of two limb cross simultaneously based on the designed two ultrasonic transmitter plethsymograph. The following additional constraints were also incorporated:

- 1. The ULP shall be capable of being attached to the crewmembers during launch and re-entry.
- 2. The ULP shall be capable of interfacing with the Experiment Data Interface (EDI) as specified.
- 3. The built-in Data Recording System (DRS) shall be capable of intermittantly recording ULP data during launch and re-entry.
- 4. The ULP shall operate on battery power.
- 5. The ULP shall be capable of operation from a DC power supply through an umbilical from the EDI.
- 6. The ULP shall be capable of operating continuously for 30 hours.
  One sample per second should be taken during this period.
- 7. The built-in DRS shall be capable of:
  - (a) Recording chord length sets determined at rates between one and ten times per second.
  - (b) Storing 1 to 4 megabits of data.
  - (c) Recording quiet period indications.
  - (d) Recording the 12 bit chord lengths generated by the ULP and the time associated with that data.
  - (e) Remotely turning on, taking data and then turning off.

## General Design Requirements:

The total package for STS ULP shall be a portable, battery powered unit with dimensions not to exceed 2.5" x 5" x 8" including the batteries. The ULP shall not exceed .8 Kgs weight. The outputs of the ULP shall interface with the EDI. An on-off switch with an indicator light shall be provided on the ULP.

# Data outputs from the ULP include:

- 1. An output of chord length D<sub>1</sub> digitally coded. This output interfaces the EDI.
- 2. An output of chord length  $\mathbf{D}_2$  digitally coded. Interfaces EDI and DRS.
- 3. Outputs of chord lengths D<sub>1</sub>' and D<sub>2</sub>' representing second leg cross-section, each digitally coded. These outputs interface the EDI and DRS.
- 4. Several outputs which can be accessed allowing internal parameters to be monitored. These outputs are defined under hardware test points (additional requirements).
- 5. Timing output and quiet period indications. These interface the EDI.

The ULP is made up of several subsystems which each perform a specified function. The subsystems should follow the specifications of the two transmitter ULVS developed at University of Kentucky and use minimization techniques to reduce size. Exceptions to these specifications are noted in the following:

## Transducers

The ULP shall incorporate six piezoelectric crystals which are used as ultrasonic transducers. These transducers will be placed across the lower limb of the subject in two sets of three. Two of the transducers in each set act as transmitters of ultrasonic pulses, and the third transducer of each set acts as a receiver. The resonant frequency of these transducers is about 2.25 MHz.

# Pinger Circuit

The pinger circuit shall provide the capability to excite a preselected ultrasonic transmitter by inducing a 100 to 150 volt, 10 nanosecond positive pulse. Modifications shall be made to eliminate the need for a -15V power supply.

# Controller Logic, Transmitter Select and Counter

The function of this subsystem is to control the flow of information throughout the system, to select a transmitter for chord length measurement and to provide a count proportional to the measured ultrasonic transit time. The basic counter frequency shall be 18 MHz for this sybsystem. This system supplies the following functions:

- 1. Provides interfacing to processor.
- 2. Provides mode selection for transducer placement and ultrasonic signal visualization.
- 3. Starts and stops counter.
- 4. Inhibits receiver section for an adjustable period of time to reduce the possibility of false output.
- 5. Selects the transmitter that is to be activated.

### Receiver Circuit and Comparator

The ultrasonic signal receiving amplifier shall have a gain of 60 dB (min) and have standard wide band video capabilities. The receiver shall be transformer coupled to the ultrasonic transducer to provide high common mode rejection. The signal to noise ratio (SNR) shall be at least 40 dB. The comparator shall be a level detector at 0.2  $\pm$  0.02 V measured at the output of the receiver amplifier.

### Processor and Memory

This subsystem provides a timing signal for operation of the ULP.

The timing signal shall be derived from a 18.0 MHz clock which will also be used for the binary count. A microprocessor, equivalent to

the 8080, shall be used for all command and store operations. The computational constants and instructions should be stored on permanent memory. Data output shall be interfaced through the processor. This subsystem shall drive a display on the ULP. The processor shall also provide command operations to the DRS.

# Additional Requirements:

#### Hardware Tests Points

A connector shall be accessible which allows the following phenomena to be monitored:

- 1. The analog signals of the transmitters.
- 2. The analog signals of the receivers.
- 3. The comparator output. (digital)
- 4. The calculated cross sectional area measurement.
- 5. The 18 MHz time base. (digital)
- 6. The Pinger voltage

#### Self-Test Button

This button should be available to activate programmed inputs into the processor for calibration and check-out.

# Quiet Period Indicator

An indicator for quiescent periods and sleep period shall be incorporated into the design of the ULP. These indications are to be imposed on the data.

### Physical Configuration

Throughout the design phase of the program, weight, volume and power consumption minimization shall be a prime consideration within the measurement, operational and environmental constraints imposed herein. The ULP electronics shall be configured in one package. An

on-off switch shall be supplied on the electronics package. This package shall be two output connectors. One output shall provide chord transit time data, timing outputs and quiet period indications. The other output connector shall be able to provide all internal information which is required herein. The package shall be easily attached to the crewmembers and be comfortable to wear.

## Range

The ULP shall have the capability to demonstrate the following characteristics:

- 1. Provide reliable measurements on calf lengths up to 18 cms corresponding to limb circumference of approximately 57 cms.
- Duration of a single cross sectional area measurement shall not exceed 50 mseconds.
- 3. The chord length measurement must be accurate to within  $\pm$  0.6 mm for a nominal chord length of 14 cms.
- 4. The digital counter shall be 12 bits.

#### Calibration

The capability to inject appropriate calibration and test signals into the STS ultrasonic plethysmograph just prior to initiation of a data gathering sequence should be provided. This will be initiated whenever the self-test button is depressed.

#### Display

A four digit display shall be provided on the ULP which displays instantaneously to the crewmembers the measured chord length of the limb cross-sectional area. The display shall be visible to the crewmembers at all time during operation.

### Battery Monitor

The ULP system shall have a battery check system and a warning if the battery power is low.

# Optimum Signal Locator

A method of detecting transducer placement for the optimum signal without using an oscilloscope shall be incorporated. This method can be an envelope detector. A blinking display can be used to show that the optimum signal is not being received.

# Ground Support Equipment

The Ground Support Equipment (GSE) to be developed in support of this program shall include as a minimum the following:

- 1. Playback equipment to allow playback of tapes recorded by the DRS.
- 2. ULP Data Interface Simulator (DIS). The DIS will be designed to emulate the ULP Data Interface with the EDI in amplitude, impedance, and format content to allow development and verification of hardware and software in, and down-stream of, the EDI. DIS electrical connectors need not be flight type.

The tape playback equipment should be supplied at the time of delivery of the first ULP prototype system containing a built-in DRS.

Delivery of the DIS will be required at least two months prior to delivery of the ULP prototype in order to allow timely development of the EDI and associated hardware and software.

Denver Research Institute, DRI, fabricated a prototype test unit for our evaluations in accordance with a test plan developed by us. Results of these tests are detailed in the "Appendix B".

# TRANSDUCER DEVELOPMENT

The use of an ultrasonic micrometer in limb plethysmography has been described in this report (page 6). This technique, ultrasonic pulses are created through shock excitation of a piezo-electric transducer. After transit through the coupling medium, a second piezo-electric transducer converts the acoustic energy into an electrical signal. If the sound velocity in the medium is known and constant, the propagation time between the two crystals is an accurate measure of the separation distance. Arrival of the pulse at receiver crystal is detected by the triggering of a comparator just above the noise level of the receiver amplifier. In the usual situation the received signal level is very much larger than the noise level; changes in comparator firing time can then be confidently interpreted as changes in tissue dimension. If pulse amplitude is not much higher than the comparator level, however, small changes signal amplitude or pulse shape will cause changes in comparator firing time. This will be falsely interpreted as a change in tissue dimension. Although these considerations are important in all applications of sonomicrometry, they are especially critical in plethysmographic applications, where very small dimension changes are measured.

To routinely obtain good signals, it is essential that well-designed transducers be available. The most critical requirement is that the beam pattern of the transducers must be sufficiently wide. Adequate signal amplitude can then be obtained with minimum care in placement of the transducers. In addition, small rotations of the transducers, which may occur during a limb expansion, will not affect

pulse shape substantially.

. . . . . . . .

Two basic types of transducers for ultrasonic dimension measurements have been described in the literature, chiefly for use in cardiac dimension measurements. One type is constructed from a disk cut from flat stock piezo-electric ceramic. Size of the disk may vary from 2 mm to 5 mm in diameter, with the smaller diameters having the wider beamwidth. A convex plastic lens is added to diverge the beam [14,15,16]. This type of transducer has the advantage that it and simple to construct. inexpensive The beamwidth transducers, however, is quite narrow, even with the addition of a Thus, careful orientation of the transducers is diverging lens. required. Small rotations of the transducers may produce relatively large changes in pulse amplitude and shape even after careful placement of the transducers on the tissue.

A second type of transducer has been more recently described [27]. This transducer uses a hemispherically shaped piezo-electric element which produces a broad beam pattern. A quarter wavelength matching layer between the ceramic and the tissue was used to maximize signal amplitude. With this approach, much less care in placement of the transducers is required. The hemispherical transducers are, however, quite expensive, and construction of the matching layer is somewhat difficult.

Although the problems encountered are similar, the transducers used in cardiometry are not ideal for use in plethysmography. In particular, transducers used for plethysmography require a rugged housing and a convenient means of attaching the housing to the skin. In this report we describe the design and construction of transducers specifically designed for plethysmographic applications, and provide


details on the performance of these units.

# Design Considerations

For attaching the ultrasonic transducers to the skin, it is convenient to use a housing of the same size and shape as that used with exercise EKG electrodes. The transducers may then be very simply attached with double-sided adhesive collars. With the size of the external housing determined, it is necessary to consider the optimum design of the transducer mounting within the housing. Both the hemispherical and the flat plate transducer designs have been mounted within the housing. We will first consider the flat transducer.

Small flat-plate ultrasonic transducers produce a wider beam, and the electrical impedance of ceramic transducers is reasonably low even for small transducers. Thus, for this application flat plate transducers should be as small as possible. For measurements across the calf, a transducer frequency of 2 MHz offers a convenient compromise between signal attenuation, which increases with frequency, and transducer thickness, which decreases with frequency. At 2 MHz the ceramic material is approximately 1 mm thick; the smallest disc which can be reproducibly cut from this material is about 3 mm in diameter. These considerations determine the size of the transducer element.

Because addition of a convex lens will diverge the beam pattern somewhat, it is desirable to include a lens in the design. If the lens is allowed to protrude into the tissue slightly, it will also serve to provide good acoustic coupling between the transducer and the tissue. The implementation of this concept is shown in Fig.21A, where a 3 mm diameter transducer is shown mounted behind a lens which






FIGURE 21. Design of Flat and Hemispherical Transducers

a 3 mm diameter transducer is shown mounted behind a lens which protrudes 3 mm into the tissue. The lens is mounted within the EKG housing to provide a space for attachment of the adhesive collars.

Hemispherical transducers may be mounted similarly within a housing of the same external shape as that used for the flat transducers. Because the hemisphere radiates in all directions, no lens is necessary; however, it is convenient to allow the transducer to protrude into the tissue to provide good acoustic coupling. Although hemispherical transducers are available in several diameters, there is no reason to expect significant differences in performance for different diameters. A more critical question is the design of the coupling layer between the crystal and the tissue.

Signal amplitude produced by an ultrasonic transducer may increased somewhat by placing a quarter wavelength thick matching layer between the ceramic and the tissue [27]. This approach routinely used in design of transducers for imaging applications [28]. In dimension gauges, however, signal levels are normally very high compared with those encountered in imaging; the wide beam pattern produced by hemispherical transducers enable strong signals to be obtained even when the crystals are badly misaligned. Thus, increasing the signal strength by several dB is of marginal value. The quarter wave length matching layer is very thin (<0.3 mm) and, thus, difficult to fabricate and is subject to damage from abrasion. The capacitive coupling between the outer conductive coating of crystal and the subject is large when a thin dielectric layer is used. The coupling can introduce high frequency noise into the system. these reasons, the developed transducers use a coupling layer with a Because this layer is more than one minimum thickness of 1 mm.

wavelength thick, the second half-cycle of the received signal is not affected by reflections within the matching layer. Yet the layer is sufficiently thin so that attentuation within the layer is negligible. The implementation of this concept with a hemisphere of 12 mm diameter is shown in Fig. 21B.

### Transducer Construction

Construction of the hemispherical receiver transducers is illustrated in Fig. 22. Three major elements are used to assemble the The plastic lens is molded from epoxy resin finished transducer. (STYCAST LN 78058, Emerson and Cuming, Canton, Mass.) in a machined After casting and curing, the back of this lens teflon mold. machined with a ball mill to fit the hemispherical transducer. transducer element is obtained commercially (Channel Industries, Santa Barbara, California) as a complete hemisphere. Teflon-insulated stainless steel leads (33 gauge, Cooner, Chatsworth, California) are carefully soldered to the inner and outer conductive coating of crystal using a zinc chloride flux and ordinary solder. The crystal is then cemented into the plastic lens. A silastic tube is over the leads. This assembly is mounted inside the pre-cast plastic exercise EKG electrode housing (Blackwell Plastics, Houston, Texas). The space between the plastic lens and the housing is sealed with noncorrosive silastic sealant (3140 RTV, Dow Corning Corp., Midland Mich).

Hemispherical transducers are available in several sizes, with the smallest being 6 mm in diameter with a resonant frequency of 3.8 MHz. Crystals resonant 2 MHz are available with a diameter of 12 mm and 9 mm. Any of these elements can be mounted within the lens

# Reciever:




FIGURE 22. Construction of Receiver and Transmitter Elements


#### structure.

The housing of the flat transmitter transducer is accomplished in a similar manner. A 3 mm diameter disc is cut from a plate of LTZ-2 transducer material resonant at 2 MHz (Transducer Products, Inc., Gophen, Conn.). Leads are attached as before and the transducer is cemented into a cavity machined into the cast plastic lens. The assembly is then cemented into the plastic housing as before.


## Transducer Performance

The external appearance of the developed transducers is depicted in Figure 23A. Fig. 23B shows the appearance of the transducers after attachment of the calf. Transducers are attached at two sites to allow simultaneous measurement of two chords. The adhesive collars provide a firm attachment which can be maintained for several days. Use of a very small quantity of acoustic coupling gel provides improved signal transmission in short-term attachment. For long-term experiments, gel is not necessary because normal skin secretions provide adequate coupling after several minutes attachment.

The most important characteristic of the transducer is the extent to which rotation off-axis reduces signal amplitude. Measured performance of the two transducer types is shown in Fig. 24. Two transducers were suspended in an acoustic test range, separated by 10 cm. In all tests a 3 mm diameter transducer was excited using a standard pulser to produce a constant acoustic signal. The transducer under test was connected to an amplifier having an input impedance of 50 ohms. The tested transducer was rotated in 5 degree increments and beam strength was determined by adjusting a calibrated attenuator to return the amplitude of the second half-cycle to the level obtained at




A. External Appearance



B. Attached to Mid-calf Region

FIGURE 23. The Developed Ultrasonic Transducers







FIGURE 24. Off Axis Performance of Transducers

O degrees. As expected the flat transducer is quite directional. The response of the hemispherical transducer is unaffected by rotation over the range of +40 to -40 degrees. Fall off of the response beyond this angle is probably due to interference within the plastic lens. The beam width of the hemisphere could be improved by protruding the lens further, but this would cause unnecessary distortion of the tissue and possible discomfort to the subject.


The directional characteristics of the transducers is further demonstrated in Fig. 25, where recorded pulses are shown. With the flat transducer the pulse shape remains unaltered with rotation, but the amplitude decreases sharply. The pulse shape of the hemispherical transducer is more complex, and changes with rotation; however, the amplitude and shape of the first cycle of the received pulse are almost unaffected by rotation. Because pulse arrival is determined by triggering a comparator on the second half-cycle of the received pulse, the remainder of the pulse is unimportant for plethysmographic applications.

In design of associated electrical circuitry, knowledge of the electrical behavior of the mounted transducers is necessary. Impedances of the two transducer types were measured in air, using an impedance meter (HP model 4815A). The flat transducer, because of its smaller size, presents a relatively large impedance (approximately 800 ohms at resonance) and behaves as a much more reactive (capacitive) load. The hemispherical transducer displays an impedance of only 80 ohms and behaves as a nearly resistive load at resonance. At very high frequencies the impedance of both transducers decreases sharply due to the parallel capacitance of the crystal.


# FLAT TRANSDUCER




10°




3/8 INCH DIAMETER HEMISPHERICAL TRANSDUCER





20°



40°


FIGURE 25. Directional Characteristics of Ultrasonic Transducers

### Use of the Transducers


Although these transducers may be used interchangeably receiver and transmitter, when the two are used in combination it is often convenient to use the flat transducer as a transmitter and the hemisphere as a receiver. Because of the high voltage and current required, changing impedance and voltage levels in the pulser requires extensive redesign of the circuit. The low impedance of hemispherical transducer at high frequencies may load the pulser some dimension gauges excessively. The flat transducer, however, will produce negligible loading. On the other hand, when the hemisphere is used as a receiver, the low impedance of the hemisphere can easily be matched to the higher impedance of the RF amplifier by use step-up transformer, producing a higher voltage input. impedance of the hemisphere allows use of the relatively long cable without significant capacitive loading of the transducer.

The use of two hemispheres provides most freedom in placement of transducers when a sufficiently low-impedance pulser is used. This is illustrated in Fig. 26A, where signals passed though a calf of 38 cm girth are reproduced. To provide a reference for location of the transduce, r the zero point of a cloth tape measure was attached to the leg at the sharp point of the tibia with the numbers increasing toward the inside of the leg. The tape was then wrapped loosely about the calf, and attached with adhesive tape. Positions of transducers were determined by recording the tape number nearest the transducer. Positions were then computed as angular deflections as if the cross-section of the calf were circular.


A hemispherical transducer was attached to the calf at 76 degrees




76/208




76/208




7.6/246



76/237



76/293



76/246

HEMISPHERE/HEMISPHERE

HEMISPHERE/FLAT

FIGURE 26. Comparison of Performance of Different Transducer Combinations

the lower edge of the tibia by palpation and then moving several cm down the calf. This location is optimum for the plethysmographic measurement: Placing the transducer further down the calf forces passage of the signal through the fibula; placing the transducer nearer to the tibia causes transducer movement to be restricted by the bone. The transmitting transducer was then placed at selected locations, as indicated in the Figure 26. Excellent signal quality was obtained as one of the transducers was moved over an 85 degree arc on the calf. The receiver amplifier is saturated on the second half-cycle of the received signal in each photograph.

In most situations the combination of a hemisphere as receiver with a flat crystal as transmitter provides very good results, although the range from which adequate signals can be obtained is smaller (Fig. 26B). Performance of this combination is due to the anatomy of the calf and the attachment procedure used. The receiver crystal (hemisphere) is first attached; the transmitter (flat) crystal is then placed against the opposite side of the leg and adjusted until maximum signal level is obtained. Thus, the directional transducer is adjusted until it is aimed at the receiver. This procedure is far simpler than would be the simultaneous adjustment of two directional transducers.

When flat transducers are used as both receiver and transmitter, signal quality is found to vary dramatically with small changes in transducer location. Small rotations of the transducer will also produce significant changes in the shape of the received pulse. These pulse shape changes are due to the narrow beam width of the flat transducers.

Both the combinations illustrated in Fig. 26 have been tested extensively both in short term physiological experiments and in long term attachment tests. Good signal quality is routinely obtained with minimal training of operators. In long-term attachment tests of up to 7 days no discomfort has been noted; signal quality is maintained so long as the transducers are firmly attached.

# Summary and Conclusions:

The development of ultrasonic dimension measurement principles into a laboratory instrument providing plethysmographic data has been detailed in this report. The developed instrument has been compared with the Whitney strain gauge and on impedance measuring device and shown to yield consistent and reliable results.

In human subjects undergoing provocative maneuvers such as venous occlusion and tilt experiments the results have been shown to be repeatable and reliable. The major aim of this project was to develop instrumentation which freed the investigator from the transducer related problems of the capacitance plethysmograph and this was achieved. Since the developed plethysmograph uses transducers on opposite sides of the limbs, the intregrity of movement is not compromised and also individualized transducer configurations are not required.

The transformation of the above instrument into a flight certified unit was carried out by Denver Research Institute. This unit has been fabricated and tested in accordance with the procedures developed in the report and is a part of NASA supplied LSLE item.

Theses, Dissertations and Publications:

The following theses, dissertations and publications resulted from the support derived from this research contract.

- i. Murali P. Kadaba "Tissue characterization using ultrasonic transmission and scattering parameters." Ph.D., 1978, Co-director - Dr. J. F. Lafferty.
- ii. An ultrasonic plethysmograph, P. K. Bhagat, J. F. Lafferty and M. P. Kadaba, Aerospace Medical Association 1979 Meeting, Washington, D. C., May, 1979.
- iii. A microprocessor based ultrasonic data acquisition and analysis system, P. K. Bhagat, M. P. Kadaba and V. Wu, AAMI 14th Annual Meeting, Las Vegas, May, 1979.
  - iv. Utility of acoustic parameters for remote characterization of tissues, P. K. Bhagat and M. P. Kadaba, Fourth Midwest Biomedical Engineering Conference and Workshop, Columbus, Ohio, May 4-5, 1979.
    - v. A microprocessor based system for ultrasonic tissue characterization, P. K. Bhagat, M. P. Kadaba, V. N. Gupta, and V. Wu, Medical Instrumentation, 14(4), 220-224, 1980.
  - vi. Attenuation and backscattering of ultrasound in freshly excised animal tissues, M. P. Kadaba, P. K. Bhagat, and V. Wu, I.E.E.E. Transactions on Biomedical Engineering, BME-27(2), 76-83, 1980.
- vii. An ultrasonic plethysmograph for space flight applications, P. K. Bhagat, J. F. Lafferty, D. Bowman, and M. P. Kadaba, Aviation. Space and Environmental Medicine, 51(2), 185-188, 1980.
- viii. Ultrasonic characterization of aging in skin tissue, P. K. Bhagat, W. Kerrick, and R. W. Ware, <u>Ultrasound in Medicine and Biology</u>, 6, 369-375, 1980.
  - ix. Estimating ultrasound propagation velocity in tissue using unwrapped phase spectra, V. N. Gupta, P. K. Bhagat, and A. M. Fried, <u>Ultrasonic Imaging</u>, 2, 223-231, 1980.
    - x. Significance of phase of an ultrasonic pulse in tissue characterization, V. N. Gupta, P. K. Bhagat and A. M. Fried, Fifth Midwest Biomedical Engineering Conference and Workshop, Columbus, Ohio, May 30-31, 1980.

- xi. A two transmitter ultrasonic plethysmograph, P. K. Bhagat, V. C. Wu, W. T. Nickell and P. K. Kadaba. Fifth Midwest Biomedical Engineering Conference and Workshop, Columbus, Ohio, May 30-31, 1980.
- xii. Ultrasonic pulse broadening due to tissue inhomogeneities, V. N. Gupta, P. K. Bhagat, A. M. Fried, and M. P. Kadaba. Fifth International Symposium on Ultrasonic Imaging and Tissue Characterization, Gaithersburg, Maryland, June 1-4, 1980.
- xiii. Vic C. Wu
  "A microprocessor based ultrasonic limb volume measurement system." MSEE, 1980, Co-adivsor with Dr. T. L. Henderson.
  - xiv. Vaikunth N. Gupta
    "Ultrasonic tissue characterization via homomorphic transform." Ph.D., 1980.
    - xv. Ruth A. G. Dyer "Classification and characterizatio of tissue pathology through ultrasonic signal analysis." Ph.D., 1980.
  - xvi. A comparison of limb plethysmograph systems proposed for use on the space shuttle, B. M. Levitan, J. Zieglschmid, L. D. Montgomery and P. K. Bhagat, Aerospace Medical Association 1981 Meeting, San Antonio, Texas, May, 1981.
- xvii. Ultrasonic properties of normal and functionally impaired renal tissue, V. N. Gupta and P. K. Bhagat. Sixth International Symposium on Ultrasonic Imaging and Tissue Characterization, Gaithersburg, Maryland, June 1-4, 1981.
- xviii. Application of orthogonal transforms for classification of myocardial backscattered signals, R. A. G. Dyer, S. A. Dyer, W. O'Connor, and P. K. Bhagat. Sixth International Symposium on Ultrasonic Imaging and Tissue Characterization, Gaithersburg, Maryland, June 1-4, 1981.
  - xix. Mark Shafer
    "The application of homomorphic processing to ultrasonic signals." MSEE, 1981, Co-director, Dr. B. J. Leon.
    - xx. Clinical applications of ultrasonic plethysmography, P. K. Bhagat, V. C. Wu, and W. T. Nickell, 17th annual meeting AAMI, San Francisco, California, May, 1982.
  - xxi. Homomorphic processing of backscattered ultrasonic signals from normal and functionally impaired renal tissue, P. K. Bhagat, M. E. Shafer, W. T. Nickell, and V. N. Gupta, Seventh International Symposium on Ultrasonic Imaging and Tissue Characterization, Gaithersburg, Maryland, June, 1982.
  - xxii. Ultrasonic quantification of skin and muscle hemodynamics, P. K. Bhagat, W. T. Nickell, and V. C. Wu, First Southern Biomedical Engineering Conference, Shreveport, Louisiana,

- June, 1982.
- xxiii. Ultrasonic characterization of acute renal failure, V. N. Gupta, P. K. Bhagat, C. E. Ott, and A. M. Fried, <u>Ultrasound</u> in <u>Medicine and Biology</u>, 8, 249-261, 1982.
  - xxiv. A practical ultrasonic plethysmograph, V. C. Wu, W. T. Nickell, and P. K. Bhagat, <u>Aviation</u>, <u>Space</u> and <u>Environmental Medicine</u>, 53(4), 375-388, 1982.
    - XXV. Transducers for ultrasonic plethysmography, W. T. Nickell, V. C. Wu and P. K. Bhagat, <u>Aviation</u>, <u>Space and Environmental Medicine</u>, 54(5), 458-463, 1983.
  - xxvi. Ultrasonic quantification of skin and muscle hemodynamics. P. K. Bhagat, W. T. Nickell and V. C. Wu, Proceedings First Southern Biomedical Engineering Conference, Shreveport, Louisiana, June, 1982.
- xxvii. A comparison of limb plethysmograph systems proposed for use on the space shuttle, B. M. Levitan, L. D. Montgomery, P. K. Bhagat and J. F. Zieglschmid, <u>Aviation</u>, <u>Space</u> and <u>Environmental Medicine</u>, 54(1), 6-10, 1983.

#### APPENDIX A

### Background Data on Transducers

# A. Piezoelectric material properties

Ferroelectric crystals are among the ceramic materials which have identical centers of positive and negative charges. Α necessary but not sufficient condition for a solid to be ferroelectric the absence of a center of symmetry of electrical charges. type of material such as quartz, either elongates or contracts upon the application of an electric field. This dimensional change is due to the effect of the applied electric field on the dipole length. Such behavior in the ferroelectric material provides a means converting electrical energy into mechanical engery. The vibrates with the frequency of applied alternating field the voltage differential. proportion to This property of ferroelectric materials of dimensional changes upon the application of an electric field has led to the electro-mechanical transducer devices, which are commonly used for the production of high frequency sound waves.

In piezoelectric material the converse effect exists, i.e., upon the application of mechanical stress to the crystal, an electric potential is produced which is proportional to the applied stress. There are 21 classes of crystals which can be classified as ferroelectric material (24). Twenty of these materials are piezoelectric, i.e., they exhibit the converse effect too. Some commonly used piezoelectric ceramics are barium titanate (BaTiO<sub>3</sub>) and lead zirconate (PbZrO<sub>3</sub>).


Ferroelectric ceramics are spontaneously polarized material.

This polarization is dependent on the temperature at which exists. The transition from ferroelectric nonferroelectric occurs at a temperature called the "Curie point". Barium titanate (BaTiO<sub>3</sub>) has a stable cubic structure above 120°C. At room temperature it changes to a tetragonal structure, thus producing a polarization due to the movement of the titanium ion  $(Ti^{+4})$ . movement of the titanium ion  $(Ti^{+4})$  introduces an electrical charge imbalance into the unit cell. There is the probability that titanium ion would be located in one of two locations along one of the major axis (x,y, or z). Since neither location is at the center of the unit cell the center of negative and positive charges are not coincidental and an electrical dipole must exist. The movement of the titanium ion to the second location by a neighboring polarized unit cell is possible, since the energy barrier between the two possible sites along one axis is low. This is known as the cooperative alignment of unit cells, which gives rise to electrical domain in the ferroelectric material as shown in Figure 27.

Upon application of an external electrical field the ferroelectric material domain, will follow the polarization of the applied field, and dimensional changes in the material will occur due to a change in the dipole orientation in the unit cell. Therefore, if an alternating electrical field was applied, vibration of the crystal will take place and ultrasound waves will be generated.

Mason (25) noted that dimensional changes in piezoelectric crystals are dependent on the polarity of the applied electric field and the crystal orientation. He deduced that the piezoelectric properties of a crystal are anisotropic and depend on its orientation.

Piezoelectric properties of a crystal are usually described by



- a) Schematic Representation of Cooperative Alignment of Unit Cells
- b) The Domain in Ferroelectric Materials




FIGURE 27. Domain Alignment in Ferroelectric Materials

certain parameters. Since these properties are direction-dependent, indices are usually used to specify each individual parameters. Some of the important parameters which pertain to the selection of a piezoelectric crystal for a specific application are listed below:

1. Transmitting Constant  $d_{ij}$ 

Defined as the strain produced in a transducer by the application of a unit electric field, and has the unit of coulomb/newton.

2. Receiving Constant gii

Defined as the electromotive force produced in a crystal upon the application of a unit pressure input, and has the unit of  $\frac{\text{volt/meter}}{\text{newton/(meter)}^2} \ .$ 

3. Deformation Coefficient hij

Obtained by multiplying the receiving coefficient by Young's Modulus of the crystal for a specified direction. It has the unit of newton/coulomb.

Coupling Factor k<sub>ij</sub>

Provides a measure of the crystal as a energy converter.

$$h = (h_{ij}d_{ij})^{1/2}$$

5. Permittivity Constant  $\epsilon$  T,S

The dielectric constant representing the permittivity of the material, and how much current is flowing through for each volt. The units of  $\varepsilon$   $^{T,S}$  are farad/meter. For this constant, the value of  $\varepsilon$  for the crystal is different in the clamped ( $\varepsilon$ <sup>S</sup>) than the free ( $\varepsilon$ <sup>T</sup>) state (25).

6. The Mechanical Factor  $Q_{\boldsymbol{m}}$ 

Determines the frequency characteristic of the crystal. This factor is independent of the geometry of the crystal, and is a

measure of the frequency band width:

$$Q_{m} = \frac{fo}{f_2 - f_1}$$

fo = Resonance frequency

f<sub>1</sub> = frequency below resonance for a reduction
 in amplitude of 3 db.

f<sub>2</sub> = frequency above resonance for a reduction
 in amplitude of 3 db.

#### B. Ultrasonic Field

When an ultrasonic transducer is excited, waves are generated and travel in the immediate medium surrounding the transducer. A pressure distribution with maxima and mimina is created in the medium. This is called an ultrasonic field. Fig. 28 is a schematic representation of the field for a circular disk transducer. Positions and amplitudes of the extrema in the ultrasonic field depend greatly upon the ratio of the surface dimensions of the source to the wave length produced.

There are two distinctive regions in the ultrasonic field which are of importance in the design of the transducer. The two regions are: The Fresnel Region (nearfield) and Fraunhofer Region (farfield).

### 1. The Fresnel Region

The Fresnel Region is also known as the near field. Its depth for circular transducer is given by:


$$x = d^2/4 \lambda$$

where:

d: disk diameter

.λ: wave length.

x corresponds to the last axial maximum pressure which is the end of Fresnel zone and the beginning of Fraunhofer zone. For a complete mathematical treatment of this region, see Hueter and Bolt (26).



Distance from the Transducer Equation for the Hearfield:  $\frac{I}{I_0} = \sin^2\left[\frac{\pi}{\lambda}\left(\sqrt{r^2 + x^2} - x\right)\right]$ 

· FIGURE 28. The Ultrasonic Field of Circular Disk Transducer

### 2. The Fraunhofer Region

The ultrasonic field generated by a circular disk transducer had the shape of a cylinder in the nearfield region. At the beginning of the farfield region the ultrasonic beam starts diverging from a uniform cylinder shape. The angle of the beam divergence from a cylindrical shape is given by:

$$\sin\theta = 0.61(2\lambda) \tag{11}$$

And for a square transducer is given by:

$$sin\theta = \lambda/t$$

Where (d) is the disk diameter and (t) is the width of the crystal. In this region the intensity of the wave is reduced sharply due to the divergance effect, and this will affect the attenuation measurement in the far field region.

### C. Measurement System

We have developed a microprocessor based ultrasonic data acquisition and analysis system for measurement of acoustic parameters of a given sample (19,23). This system was used for measurement of transducer beam profile and central axis intensity as a function of distance.

### Instrumentation details:

The main components of the developed pulse-echo system are shown in a block diagram in Fig 29. The system consists of a pulse generator, radio frequency (RF) receiver amplifier, calibrated attenuator, and a peak detector. A range gate enables the selection of a particular segment of the reflected signals for processing. All





Figure 29. Data collection and analysis block diagram

of the above are housed in a single unit (Aerotech Labs, 1 UTA-3). input receiver sensitivity is 5 mV with a frequency bandwidth of MHz and a lower frequency cutoff at 400kHz. The receiver gain is 40dB with an output impedance of 50  $\Omega$ . A second set of pulser, receiver, range gate (Metrotek, Inc.2) was utilized in some experiments primarily due to the higher value of excitation (adjustable up to 230V) delivered by the pulser. Acquistion of a microprocessor (Altair 8800), a high speed A/D convertor (Biomation 8100) and an off line minicomputer (PDP 11/34) facilitates digitization and storage of rf Accurate transit time signals on a disk for further processing. measurements can be accomplished using a 225 MHz universal counter timer (Textronix 7B15). Spectral analysis of data can be performed in analog fashion using spectrum analyzer (Tektronix 7L13) or digitally using standard Fast Fourier Transform (FFT) routines available on the DEC LAB applications package.

# Mechanical Apparatus

The mechanical apparatus used in this project, shown in Figure 30, was designed and assembled in our laboratory. The acoustic test tank is made of clear, 1.25-cm thick plexiglass, and is lined with a sound-absorbing Styrofoam sheet to prevent any extraneous reverberations caused by multiple reflections. The test tank is filled with degassed distilled water, which acts as the coupling medium between the ultrasonic transducer and the specimen. The transducer holders are made out of plexiglass and are mounted on a sterotaxic positioning device, which provides several degrees of freedom for accurate positioning of the transducer. Thick plexiglass slabs of various geometries were

Subsidiary of KB-Aerotech, Lewistown, PA. Metrotek, Inc., Richland, WA



used as reflectors. A 1-rpm synchronous motor is used to control the angular orientation of the transducer with respect to tissue sample.

The mechanical scanner consists of two motor-driven slide assemblies capable of movement in a longitudinal direction, as well as in a transverse direction. The linear velocity of the slides may be separately controlled. The slides have a maximum displacement of 30 cm in the longitudinal direction and 5 cm in the transverse direction. Analog position indicators were incorporated to indicate the position of the scanning transducer with respect to a reference coordinate frame. The whole assembly is made of Duralmin. A calibrated hydrophone (Mediscan, Inc.) with a crystal diameter of 0.7 mm is used as receiver to scan the beam. The scanning is carried out at various axial distances from the transmitting transducer and covers the entire transmitting cross section.

# D. Piezoelectric transducer material used in this study.

The transducers used in this project were made from two types of piezoelectric material, LTZ-2 and LTZ-5. (Transducer Products, Inc.). The material is available in form of squares (1",1",t) where t is the thickness which varies to yield different resonance frequencies. The specification for the material are as follows:

| PROPERTY                                   | <u>LTZ-2</u>         | <u>LTZ-5</u> |
|--|----------------------|--------------|
| g <sub>33</sub><br>volt/m/N/m <sup>2</sup> | 25.5*10 <sup>3</sup> | 35.2         |
| d <sub>33m/volt</sub>                      | 405*1012             | 155*1012     |
| Q <sub>m</sub>                             | 75                   | 500          |

For application in ultrasonic plethysmography either material type has been observed to be suitable.

#### APPENDIX B

Test Procedure and Results Using Prototype ULP

### Introduction

The Ultrasonic Limb Plethysmograph, fabricated by Denver Research Institute, was received by the Wenner-Gren Research Laboratory (WGRL) on May 20, 1981. Tests were carried out in accordance with the test plan developed by WGRL in consultation with Clete Booher of NASA. Personnel involved in the testing program were: P. K. Bhagat, V. C. Wu, W. T. Nickell, Dave Brewer, Mark Shafer, and WGRL support personnel. Results of these tests are presented in the following pages.

#### A. PHYSICAL EXAMINATION

- I. External Dimensions (LxWxH): 18.2 cm x 11.4 cm x 5.08 cm
  Weight (kg): 1.173 kg Batteries (kg): 0.114 kg
- II. Connectors. Are all external connectors:
  - a. Sufficiently secure to ensure reliable electrical contact? Yes

العامل المحافظ المواقع المحافظ المحاف

### Comment:

A marker is needed on the female connector to facilitate correct insertion.

b. Rapidly disconnected in an emergency? X

#### Comment:

Connection appears to be tight. It should be verified by NASA.

III. Are there sharp edges on the case? X

#### Comment:

When cassette door is opened, sharp edges are exposed on the door and the locking mechanism.

IV. Check rigidity and security of indicated parts:

a. Case

O.K.

b. Batteries

X

#### Comment:

- 1. Banana jack needed tightening.
- Not easy to change batteries.
- 3. Needs a cover to prevent accidental operation of the switch on the cartridge. Recommend a supply of spare cartridges be provided with the unit.
- c. LCD Assembly

O.K.

d. Electronic Cards

X

#### Comments:

- PC board need cleaning.
- 2. Unused switch on the tape recorder needs to be covered.
- 3. Battery and tape recorder power supply wires need to be harnessed or secured.
- 4. The real time clock battery holder spring needs to be modified so as to prevent short circuit of the battery.
- V. Shake the unit. Are any unusual sounds produced? Yes

#### Comment:

Cassette is not securely mounted. Cassette strap ring should be removed.

VI. Is the LCD display readable? Yes

# VII. Tape Drive

- a. Does the tape drive motor produce an abnormal noise level? No
- b. Operate the tape removal and insertion mechanism. Is performance satisfactory? Yes
- c. Is 'rewind' switch operation satisfactory? See comments

#### Comments:

- 1. The door does not provide easy access to the rewind switch. Suggest additional clearance when door is opened (modify hinges).
- When the door is closed, movement of the recorder switch to "Stop" position should be prevented.
- 3. Operation of door closing mechanism is difficult.
- IIX. Check materials included in prototype against materials list.

# Mechanical Listing

- 1. Toggle switch movement shown in drawing ED-14843 is 90° off prototype position Change the drawing.
- 2. Prototype fitted with Phillip screws. Drawing ED-14843 shows straight screws.
- Items in the drawings should specify all parts. For example, Item 7 in ED-14843 should specify both push button switches.

Needs electrical parts listing and materials used list.

IX. Can all indicator lights be seen under normal ambient light levels? Yes Color as required? Yes

### Comment:

Ϋ.

When the tape is full, there is a warning sound. This sound should be approved by NASA for operational convenience. Audibility of this sound under Spacelab condition should be checked by NASA.

# B. ELECTRONIC TESTS

- I. Circuit diagram review. Carefully review the units circuit diagram.
  - a. Is item complete? Yes
  - b. Is circuit satisfactory? See comments

#### Comments:

- i. Symbols not consistent example LED
- ii. Core and size of transmit transformer not defined.
- iii. Receiver transformer ratio not defined. Suggest using the receiver transformer for impedance matching between transducer and amplifier input stage.
  - iv. Pin 16 of hybrid ED-15236-2 is not defined.
  - v. The signal quality test light is amber instead of red as shown in drawing.
  - vi. In EC-15033 and ED14401, a bigger DIP package than needed is used.(!)
- vii. Size and type of plug should be defined for DRS power control jack and DRS "MIC" jack (EC-14402).
- viii. Capacitor on pin 7 of power supply hybrid not labeled (EC-14402).
  - ix. Pin 17 test point not identified in the drawing EC15019.
- II. Functional Test. Place a transducer set in a water bath. Verify that instrument works, using LCD display as output.
  - a. Stability

O.K.

b. Accuracy

O.K.

#### Comments:


- a. Up to 1 mm fluctuation in display.
- b.  $24 \pm 1$  counts/2 mm
- III. Quality of output signals. Conduct these measurements with power supply voltage at lowest specified value. Photograph signals using oscilloscope and camera.
  - a. Output pulse. Measure with 50 ohm termination.

Amplitude? <u>5V.P-P</u>


Rise Time? 35 nS

| b. | Frequency content. Determine pulse.    | e power sp   | ectrum of     |
|----|--|--------------|---------------|
|    | See Attached Drawing                   | Amplitude?   |               |
|    |  | Rise Time?   |               |
| c. | Output pulse. With transducer at       | tached.      |               |
|    |  | Amplitude?   | 67V.P-P       |
|    |  | Rise Time?   | 50 ns         |
| ď. | Digital output. Measure using enation. | xpected logi | c termi-      |
|    |  | Amplitude?   | 5V.P-P        |
|    |  | Pulse Shape  | ? <u>0.K.</u> |
|    |  | Noise? O.K   | <u>.</u>      |

See next pages for graphs and pictures.



Transmit Transformer output with 500 (using Biomation)



With Transmitter (Captured on Biomation)

FIGURE 31: ULP Transducer Output Pulse Shape




FIGURE 32. Digital Output in EDI Mode (2400 Baud)






FIGURE 33. Output Pulse with 50  $\Omega$  Termination





FIGURE 34. Amplitude Spectrum of the Output Pulse

OUTPUT PULSE WITH TRANSDUCER ATTACHED



91

| IV. | Receiver and Trigger.  | Conduct t | hese | tests | with | lowest |
|-----|------------------------|-----------|------|-------|------|--------|
|     | permitted supply volta | ge.       |      |       |      |        |

a. Receiver amplifier Gain

53db

Acceptable?

b. Noise level of receiver.

50mVRms

Acceptable?

c. CMRR of Receiver

>20dB

Acceptable?

226GB

d. Threshold stability of trigger.

# Acceptable?

Yes

e. Input Impedance of receiver. Measure input impedance of receiver between 450 kHz and 10MHz.

Data enclosed. Resonance occurs between 1.5 and 2 MHz (Phase Change) - Measured using vector impedance meter.

f. Frequency response of receiver. Determine 3dB point.

50kHz and -5.5MHz

Comment:

Measured using differential sine wave input.

g. Determine range of transmitter - receiver combination, as compared with original ULVS unit.

Minimum - 6.1 cms, Maximum - 24 cms

Comparable? Yes, but the minimum is higher than the original ULVS unit.

# ULP RECEIVER IMPEDANCE

| FREQ.<br>(MHz) | IMPE<br>(Oh | DANCE   |         | ASE<br>GLE |  |
|----------------|-------------|---------|---------|------------|--|
| •              | INPUT 1     | INPUT 2 | INPUT 1 | INPUT 2    |  |
| 0.5            | 370         | 380     | +90°    | +90°       |  |
| 1.0            | 1.05K       | 1.05K   | +86°    | +86°       |  |
| 1.5            | 3.0K        | 4.3K    | +60°    | +56°       |  |
| 2.0            | 3.0K        | 2.9K    | -62°    | -65°       |  |
| 2.5            | 1.4K        | 1.35K   | -76°    | -77°       |  |
| 3.0            | 950         | 920     | -80°    | -81°       |  |
| 3.5            | 720         | 710     | -82°    | 82°        |  |
| 4.0            | 600         | 580     | -82°    | -83°       |  |
| 4.5            | 500         | 480     | -82°    | -83°       |  |
| 5.0            | 430         | 420     | -82°    | -83°       |  |
| 5.5            | 380         | 370     | -82°    | -83°       |  |
| 6.0            | 340         | 335     | -82°    | -83°       |  |
| 6.5            | 310         | 300     | -82°    | -83°       |  |
| 7.0            | 280         | 280     | -82°    | -83°       |  |
| 7.5            | 260         | 260     | -82°    | -83°       |  |
| 8.0            | 240         | 235     | -82°    | -83°       |  |
| 8.5            | 220         | 215     | -82°    | -83°       |  |
| 9.0            | 200         | 200     | -80°    | -82°       |  |
| 9.5            | 195         | 195     | -80°    | -82°       |  |
| 10.0           | 185         | 180     | -80°    | -82°       |  |

#### Comment:

- Due to the design of the testing pad additional noise is introduced into the receiver output. For this reason, the numbers listed above do not reflect actual values.
- 2. The actual receiver gain is lower than that specified in the end item specification. However, operation appears to be satisfactory.
- V. Function switches. Check function of switches.

| a. | Time Set        | O.K. |
|----|-----------------|------|
| b. | Time on/off     | O.K. |
| c. | Fast/slow       | O.K. |
| đ. | Self-test/Reset | O.K. |
| e. | Power           | O.K. |

#### Comment:

- i. Operation of 2 or 4 chord switch has no effect until 'Reset' is operated.
- ii. Access to function switches is difficult due to their location under the tape recorder. Suggest vertical mounting of the switches for ease of operation.
- VI. Automatic signal quality check.

Does automatic signal quality check system perform properly? Yes

Describe measurements made to verify operation of this function:

Transducers 100 mm apart were placed in a water bath.

- a. One of the transducers was rotated until poor signal quality was apparent. The warning light indicated this condition correctly.
- b. Keeping one of the transducers fixed, the other was given sinusoidal displacement at varying rates. When the amplitudes and frequencies of 5 mm, 0.5Hz; 2.5 mm, 0.7Hz and 1.0 mm, 1.6 Hz were exceeded the warning light was 'ON' even though the signal quality was good.

The warning light sometimes is on during physiological tests even though the quality of signal is good. Suggest changing the test range of signal quality.

| VII. | Power Consum | aption. | Measu  | re avera | ige s | steady | state  | power |
|------|--------------|---------|--------|----------|-------|--------|--------|-------|
|      | consumption  | using s | series | ammeter  | and   | volt   | meter. |       |

| a. | Slow | Mode | 179-182 | ma, | 7V | 2 | chord |
|----|------|------|---------|-----|----|---|-------|
|    |      |      | 182-186 | ma. | 7V | 4 | chord |

c. Recorder Running 167 ma, 7V 2 chord 120 ma, 3.5V 2 chord

For 4 chords, the numbers are essentially the same.

d. Battery Life Under Normal Operation from full charge. 33 hrs.

Comment:

Out of spec conditions on voltage were observed on two occasions. Suggest recheck of power supply hybrid.

# IIX. Time Base Accuracy

- a. Clock Frequency 18.0042 MHz, 60.0014 Hz
- b. Determine clock Accuracy by test against official standard.

Drift was 6 seconds for four days and 37 seconds for twenty days when tested against WWV.

|         | Initial Readings | Standard    | ULP |
|---------|------------------|-------------|-----|
| l week  |                  | <del></del> |     |
| l month |                  |             |     |

# IX. Tape Operation

a. Quality of Recorded Signal

Amplitude O.K.
Noise Level O.K.

See attached data on page

b. Interface to LSLE computer. Comment on any problems encountered:

Interference between EDI data ouput and receiver input was observed, however, rearrangement of external wire eliminated this problem. We feel that EDI output sometimes feeds signals to the receiver amplifier during data transmission.

c. Data Reliability. Record a string of at least 5,000 data points. Play recorded data into LSLE computer. Determine number of transmission errors which occur.

Slow Mode Error Rate

l error

Fast Mode Error Rate

NONE

Comments:




FIGURE 36. Tape Recorded Signal Played Back with 1K Load and Capture by Biomation A/D Converter

# C. USE TESTS

| I. | Ease of use.  | Attach the  | unit to an   | experimental   | subject. |
|----|---------------|-------------|--------------|----------------|----------|
|    | Instruct subj | ect to oper | ate controls | s, change tap  | es, con- |
|    | duct simple e | xercises. a | nd remove ar | nd reattach ti | he unit. |

a. Accessibility of Controls

O.K.

b. Visibility of Indicator Lights

O.K.

c. Ease of Changing Tape and Battery

See section A, IV b of this report.

d. Difficulty in detaching or reattaching cables and connectors

See section A, II of this report.

e. Comfort

Comments:

Need 'pouch' to comment on comfort.

- II. Provocative Tests. Conduct Venous Occlusion and tilttable tests to produce changes in limb volume.
  - Transmit data directly into LSLE microcomputer

Quality of Data

O.K.

Interpretation of Data

O.K.

#### Comments:

More noise spikes on 2nd chord. Does not seem to be on transducer. The problem is intermittent. Data enclosed.

b. Record data on tape.

Quality of Data


O.K.

Interpretation of Data

O.K.

### Comments:

Data enclosed on pages 99, 100 and 101.



:

FIGURE 37. Time Course of Dimension Changes Under 50 mm Hg Venous Occlusion Procedure





FIGURE 38. Time Course of Limb Dimension Changes Under 50° Head Up Tilt



ULP Data Played Back From the Recorded Venous Occlusion (50 mmHg) Data FIGURE 39.

III. Routine Monitoring. Use tilt-table to produce leg volume data over a 30 minute period.

. Quality of Data

O.K.

Ease of Interpretation

O.K.

Comments:

Data enclosed on page 103 .

IV. Automatic signal check. Use automatic signal check feature to place transducers on subjects leg. Verify quality of signal using oscilloscope.

# trials/# correct

O.K.

Comments:




FIGURE 40. Passive Monitoring of Data Through a Tilt Table Experiment




FIGURE 40. (Continued)

### TEST EQUIPMENT USED IN EVALUATION OF ULP

- 1. Biomation, Model 8100, Transient Recorder
- 2. Tektronix, 604, Monitor
- 3. Tektronix, SC 502, Oscilloscope
- 4. Hewlett-Packard, 7402A, Strip Chart Recorder
- 5. Fluke, D802 Multimeter
- 6. Fluke, 8020A Multimeter
- 7. ULVS, Wenner-Gren Research Laboratory
- 8. OHAUS, Harvard Trip Balance
- 9. Wavetek, HF Sweep Generator, Model 144
- 10. Tektronix, 7633 Multimode Storage Oscilloscope, with plug-in units, 7A26 Dual Time Trace Amplifier, 7D15 225 MHz Universal Counter/Timer and 7B53A Dual Time Base
- 11. Hewlett-Packard, 4815A RF Vector Impedance Meter
- 12. Tektronix, 4051 Graphics System With 4051E01 Rom Expander
- 13. Tektronix, 4662 Interactive Digital Plotter
- 14. Tektronix, P6202 Probe
- 15. Tektronix, C-5B Oscilloscope Camera
- 16. Hewlett-Packard, 5303B Frequency Counter
- 17. Velmex, Series B2500, Positioning System
- 18. Velmex, Model 311, Control Unit
- 19. La Berne Tilt Table
- 20. MTS 442 Controller
- 21. Cromemco, System Three, Computer
- 22. MTS Servoram
- 23. LAMBDA LH-124 FM, D.C. Power Supply
- 24. RCA WP-702A, Dual D.C. Power Supply

# Conclusions and Recommendations

The basic design of the ULP unit is sound. Most problems encountered in our testing program were operational difficulties which can be fixed by either modifications of software and/or development of proper procedures. The following specific recommendations are made:

- 1. RF signal cannot be seen on a conventional oscilloscope because of the low repetition rate. This makes transducer placement very difficult without the original ULVS. Software fix is needed to provide at least 100 sweeps/sec.
- 2. We have had some problems with the automatic signal quality check. It is sometimes 'ON' when the signal is good and 'OFF' when the signal is bad. Resetting the levels might cure the problem. We suggest calculation of power requirement of the Automatic Signal check. The signal quality check could be made only once every 10 sec., for example if this would save power. In addition, the maximum sampling rate could be increased.
- 3. Polarity of the leads is a problem. The first 1/2 cycle of a "good" signal is normally small. It is desirable to trigger on the the 2<sup>nd</sup> half cycle. With the ULVS this problem is handling by reversing polarity. The leads should be designed so that the polarity can be set using an oscilloscope and recovered when leads are reattached.
- 4. The receiver transformer needs further work. The turns ratio should be set to match transducer impedance to receiver impedance. The transformer does not behave in an ideal manner over the important frequency range. There is a resonance between 1.5 and 2.0 MHz. We suggest consultation with WGRL or purchase of

- commercial units. Improvement of the transformer design will probably improve the receiver gain significantly.
- 5. Transmitter pulse width is long compared with the natural period of the transducer; the pulse is bipolar. This might cause problems if a higher frequency transducer were to be used.
- 6. The size and weight of the unit are slightly above specifications.

  We offer the following suggestions for reducing these values:
  - a. The present case is made of very thick material. Size and weight could be reduced significantly by using a less substantial case.
  - b. The tape recorder case, the replay electronics, and the speaker are superfluous. Space and weight could be saved by removing these components.
  - c. A more compact shape might be produced by careful arrangements of parts. Unfortunately, the length and width of the unit are set by the size of the PC board. Placement of the circuit on two PC boards would allow a more compact shape to be produced.
- 7. The gain of the receive amplifiers is slightly lower than that of the original ULVS. Increasing the gain by a factor of two would improve device performance.
- 8. Turning the unit 'OFF' causes loss of data back to the last tape recorder dump. This will be an operational inconvenience in some cases. The software should be written so that the ULP stops recording data when the transducers are disconnected. This would allow the power switch to remain 'ON' for long periods without draining the batteries significantly.

#### APPENDIX C

Description of Individual Circuit Blocks of Two Transmitter ULVS

A. Processor, Memory and Display


block (Fig. 41) consists of one Intel 8080 CPU, two EPROMs for storage of instruction codes and constants, two 2111 RAM storage of generated data, one 8228 data bus and status controller, and an 8224 clock generator which provides the working clock at 2 MHz for the CPU and also supplies the clock frequency for It also contains a 6820 Motorola PIA which the counter at 18 MHz. drives four BCD to seven segment decoders/drivers 7448, which in turn drive 4 DL-704 common cathode LED's where the results are displayed. The 5th digit of the LED display is wired so that it will show blank "1" only; this digit is driven by an npn-transistor which controlled by line CB2 of the 6820 PIA. Two 7400 NAND and one NOR gates are used to provide the required address decoding. Thev also provide suitable enabling signals. Because the 8080 has memory reference and I/O reference instructions, the circuit arranged so that the 6820 occupies I/O addresses 04 to 07 of the (address 00 to 03 is reserved for the use of the other 6820 which will be discussed later). The two 2708 ROM chips use memory addresses from 000 to 7FF Hex. Two 2111 RAM chips then occupy memory addresses from Hex. The decoding and enabling logic of each Hex to 8FF peripheral chip is implemented as follows:

Chip Select of 2708 (1) = All + MEMR . Al0

Chip Select of 2708 (2) = All + MEMR . Al0

Chip disable of 2111 (1), 2111 (2) = Al1

Output disable of 2111 (1), 2111 (2) = MEMR



Schematic of processor and memory for the two transmitter ULVS FIGURE 41-

June 27

Read/Write of 2111 (1), 2111 (2) = MEMW

Chip Select 0 and 1 (CSO and CS1) of 6820's = 1

Chip Select 2 (CS2) of 6820 (2) = A2

Chip Select 2 (CS2) of 6801 (1) =  $\overline{A2}$ 

Enable of 6820's = IOW + IOR

Read/Write of 6820's = IOR

Where AO to All are address lines from the CPU, and MEMR, MEMW, IOR, and IOW are signals from the 8228. These signals represent memory-read, memory-write, input/output read, and input/output write respectively. Note that in order to simplify the decoding circuit, non-fully-decoded addressing is used.

The Intel 8080 microprocessor serves as the CPU of the system. The two clocks,  $\beta_1$  and  $\beta_2$ , the reset, ready, and synchronizing signals are all provided by the 8224 clock, which controls the timing for the entire system. An 18 MHz crystal is connected at pins 14 and 15 of the clock and the resulting oscillatory output is available at pin 12 to run the counter in the next block (see Figure 42). An external reset button is connected between the Reset Input (Resin), Ready Input Rdyin), and ground to allow manual reset of the system. The status strobe (STSTB) is connected to the 8228 Data Bus (Pin 1, STSTB) to synchronize its operation.

The 8228 Data Bus and Status Controller is used to increase the driving capacity of the CPU. It generates all signals required to directly interface RAM, ROM, and I/O components. It also provides isolation of the 8080 data bus from memory and I/O which allows for the optimization of control signals. The isolation of the bus driver also provides for enhanced system noise immunity. Data lines D0 to D7 of the 8080 are tied to D0 to D7 respectively of the 8228. The Write




FIGURE 42 - Schematic of controller logic, transmitter select and counter for two transmitter ULVS

(WR), Data Bus in (DBIN) and the Hold Acknowledge (HLDA) lines of the 8080 are tied to the corresponding lines, WR, DBIN, and HLDA of the 8228. The DBIN signal indicates to external circuits that the data bus is in the input mode. The HLDA appears in response to the Hold signal and indicates that the data and address bus will go to the high impedance state. Line 22, BUSEN, is the bus enable input. It is tied to ground to keep the data bus always enabled. The output of the 8228 are lines DBO to DB7, the extension of 8080's data lines. These data lines are then available to the peripheral components. This provides a two-way data path with high fan-out between CPU and external components. Four other output lines of the 8228, Memory Read (MEMR), Memory Write (MEMW), I/O Read (IOR), and I/O Write (IOW) are available to control read/write operations in other components.

Memory capabilities of the system include two 256 x 4 RAM and two lk x 8 ROM chips. Input/output data lines of the two RAMs (pins 11 - 14 of both chips) are tied to the Data lines DBO-DB7 of the 8228 to allow data transfer between the CPU and RAM. Each RAM has its address line AO-A7 tied to AO-A7 of the 8080 to enable the CPU to store and read data at specific addresses in RAM. Read and Write enabling signals come from the MEMR and MEMW lines of the 8228. Input/output data lines of each ROM (pins 9-11, 13-17 of each chip) are also tied to the Data lines DBO-DB7 of the 8228 to allow data transfer between the CPU and ROM. Each ROM has its address lines AO-A9 tied to AO-A9 of the 8080 to enable the CPU to read data from specific addresses in the ROM. Each memory chip is located at the addressable location stated above.

A 6820 Peripheral Interface Adaptor (PIA) is included to

facilitate efficient data transfer of the computed limb cross-sectional area to the display. Data lines of the 8228 are tied to the corresponding data lines to the PIA and when the area has been computed, the PIA reads in the data. The data is then output from lines PAO-PA7 to the 7448 BCD-to-seven segment decoders, whose output goes to the seven segment displays. Read and write signals for the PIA are supplied by the IOR and IOW lines of the 8228 respectively. A Reset signal is supplied by by 8224 clock. The PIA is located at the addressable memory locations stated above. Switching in all circuits is controlled by the program in 8080.

## B. Controller Logic, Transmitter Select, and Counter

The function of this block (Fig. 42 ) is to control the flow of information throughout the system, to select a transmitter for chord length measurement and to provide a count proportional to measured ultrasonic transit time for further computation. It contains a 6820 PIA which drives three cascaded 4-bit 7493 binary counters, a CD4052 analog switch, a 7474 D-type flip-flop, a 74121 one shot, and 7400 NAND chip for logic level converting and gating between the binary counters, flip-flop, one shot and the analog switch.

The PIA in this block is used to interface the binary counters with the CPU, and it is connected almost in an identical manner to the previous PIA (Fig. 41). Data line DO-D7 (pins 26-33) as well as reset, enable, read/write, RSO, and RS1 (pins 34, 25, 21, 36, and 35) are all connected as described earlier. The only difference is the chip select 2 (CS2) pin. This allows the CPU to select one or the other of the PIA's. Lines PAO-PA7 and PBO-PB3 are connected to the output of the 7493 binary counters as shown in Fig. 42, in order to allow the CPU access the binary count. The three counters are

connected in cascade such that the least significant bit of the count is output at pin 12 of the 7493 #1 and the most significant bit is output at pin 11 of the 7493 #3. The counting frequency is supplied of by the 8224 clock (Fig. 41 ) and is connected to pin 14 of the 7493 #1. Counting frequency is 18 MHz, hence the accuracy of time measurements is 55 nsec. The counters are enabled by a 7474 D-type flip-flop.

The 7474 is set by the transmitter enabling pulse from the CPU through the PIA at line PB7. Pin 4 of the flip-flop is the reset. Initialization of the counters is provided by the output, Q, of the flip-flop. Q and the 18 MHz pulse from the 8224 clock are input to a 7400 NAND gate, and the output is used to enable the counter and control its counting frequency. Pin 1 of 7474 is the clear input, and it is connected to the receiver's comparator output in order to stop the counters when the signal is received.

Since the receiver amplifier has high gain, the possibility of false output from the comparator due to the presence of low level interference may exist during the activation phase of the transmitted signal. This possibility is avoided through use of a 74121 one shot which inhibits the receiver comparator output during this period. This one shot is enabled by the transmitter enabling pulse from the CPU through the 6820 at line PB7. The output of the 74121 is set at about 10 usec. This time period is of sufficient length to prevent false firing of the comparator (set as a threshold detector).

A CD4052 analog switch is used for selection of either one of the two transmitters for activation. The transmitter enabling pulse from line PB7 of the 6820 (1) is input to the 4052 where it is multiplexed


to output at either  $X_O$  (pin 12) or  $X_1$  (pin 14). Two transistors (2N222) at the output are used to convert the logic levels from the analog switch to a 15 V pulse to activate the ultrasonic transmitting circuit (pinger). S1 (pin 9) and SO (pin 10) of the 4052 are the multiplexer select lines. S1 is grounded (Logic Level "O") since only two of the outputs are used and SO is connected to the mode switch.

A two-pole three position slide switch is used as a mode selector. In the first two positions a ground level is sensed by line PB4 of the 6820 and notifies the CPU to send out approximately a 2 Khz pulse test signal to either one of the transducer enabling inputs of the transmitting "pinger", depending on the position of the switch. This in turn will provide a continuous pulsed ultrasound output so that the operator can place the receiver transducer at a proper position. In position 3, the run mode, the CPU takes over and selection of the transmitter is controlled by the stored program. Power Supply Circuits

Power for all clocks of the system is supplied by the TM515 power module. Fig. 43 shows the designed voltage regulator circuits to convert this power source into the required 5V, 12V, and 15V required by the system. The three circuits operate as follows:

Figure 43A produces +15V and +12V. A regulated voltage source provided by a 15V Zener diode and a pnp power transistor produces the -15V supply. A second similar voltage source using a npn transistor produces the +15V supply. The +15V supply is also used as the input to an LM 7812 +12V precision voltage regulator which produces the 12V supply. Capacitors are used to filter out any AC ripple.

Since the LED displays do not require a precise voltage source, a simple voltage dividing network, shown in Fig. 43B, is used to



\*

NAS 9 - 15452 Principal Investigator: P. K. Bhagat June 27, 1979

WENNER GREN

LABORATORY RESEARCH

A ALL THE RESISTOR USED HERE IS IW

43- Schematic of power supply circuits used in two transmitter ULVS

FIGURE

produce its 5V supply. Line voltage is stepped down to 25V through a transformer. A pnp power transistor and three resistors are used to produce the voltage divider which produces 5V at the output of Fig. 43 B. In order to prevent the 100 mf capacitor from discharging back into the transformer, a diode is included at the input. Capacitors are used at the input and output to filter out any AC ripple.

Two +5V precision voltage regulators (LM 7805) and one -5V precision voltage regulator (LM 7905) are used in Fig. 43B to produce the 5V supplied. Line voltage is stepped down to 15V through a center tapped transformer. The center tap is grounded. Full wave rectification of the secondary voltage across the end terminals of the top portion of Fig. 43c is accomplished by tow diodes. A regulated voltage source, provided by an npn transistor and a 9V Zener diode with input and output capacitors for filtering out ac ripple, is used input to the two 7805's. Half wave rectification of as secondary voltage across the end terminals of the lower portion of is accomplished by one reverse-biased diode to allow only negative voltages to pass. This is used as the input to 7905, which produces the -5V supply. Transistors  $Q_1$  through  $Q_2$  are power transistors supplied by the TM 515 power module. All resistors are 1 watt.

#### SOFTWARE IMPLEMENTATION

The operation of the Ultrasonic Plethysmograph is controlled by the 8080 microprocessor and the program implemented in the 8080 provides the process for determining limb cross-sectional area. To fully understand the operation of this device it is necessary to understand the software implementation. A flow graph of the program is shown is Fig. 44.

The program counter within the 8080 microprocessor chip is initialized through a SPST pushbutton switch located on the front panel. Automatic power-on-reset of the program counter is provided through an R-c charging network. The program begins initialization of the 6820 PIA. The stack pointer is then set up followed by the choice of operating mode. In the testing mode a 2 Khz pulse repetition frequency (prf) is used to excite one of transducers based on the position of the selector switch. This feature enables one to observe the received signals on oscilloscope screen and adjust the transducer locations for optimum results. Once the optimum transducer sites have been located, the selector switch is moved to position 3, run mode, indicating that limb area computations may now be carried out. It begins by resetting the binary counters upon activation of transducer 1. The program then calls subroutine "NOTYET" shown in Fig. 45. This subroutine takes the binary count, which is the time required for the ultrasonic pulse to traverse from transmitter to receiver, and converts it into a corresponding chord length, dl. This information is stored in the RAM. In a similar manner the second chord length data, d2, is obtained and stored. The limb cross-sectional area is computed and




FIGURE 44. Flow Graph of Software Implementation




FIGURE 45. Subroutine "NOTYET"

displayed on the front panel LED readout. After a delay (firmware adjustable, normally set a 1 sec), the program returns to the testing mode and a second sequence of events begins, unless the selector switch is moved to positions 1 or 2.

A complete listing of software is included in the following pages.

## SOFTWARE PROGRAM LISTING

ULTRASONIC LOWER LINB AREA COMPUTATION. \* THIS IS A PROGRAM TO COMPUTE AREA OF A SEGMENT OF THE LOVER LIMB UNDER MEASUREMENT. COMPUATTIONS ARE DONE BY UTILIZING FLOATING POINT ROUTINES WHICH ARE ADOPTED AND ; HODIFIED FROM THE INITIAL VERSION OF INSAI ; 8K BASIC SOURCE LISTING. THESE ROUTINES ARE : LISTED FOR REFERENCE. PARAMETERS ARE DEFINED AS FOLLOWS: : D1= DISTANCE BETWEEN 1ST TRANSMITING TRANSDUCER AND THE RECEIVING TRANSDUCER. D2= DISTANCE BETWEEN 2ND TRANSMITING TRANSDUCER AND THE RECEIVING TRANSDUCER. D= DISTANCE BETWEEN TWO TRANSHITTING TRANSDUCERS. WHICH IS PRESETED AS 2.5( CM). PROGRAM IS STORED IN TWO 2708 EPROMS AND STRAT AT LOCATION O. GENERATED AND TEMPERARY DATAS ARE STORED IN TWO 4 BY 256 RAM WHICH ARE LOCATED AT 800H. FLOATING POINT DATAS ARE STORED IN 4 CONSEQUTIVE BYTES WITH THE 1ST BIT OF THE 1ST BYTE

OF THE 1ST BYTE IS EXPONENT AND THE FOLLOWED

: BE THE SIGN OF THE MANTISA AND 2ND BIT OF THE 1ST

ORG

3 BYTES ARE MANTISA.

; set the stack pointer at 43000(980H).
; Locations 800H to 97FH are saved to

STORE GENERATED AND TEMPERARY DATAS.

BYTE BE THE SIGN OF EXPONENT. THE REST

LXI SP,4300Q

; initialize plas(notorola 6820).
; pia(1) Channel a is set as input.
; pia(1) first 2 bits of channel b
; is set as output others as input.

```
; INPUT ARE FROM BINARY COUNTER, OUT
 PUT ARE TO THE CONTROLING LOGIC.
  PIA(2) CHANNEL A AND B ARE SET AS
; OUTPUT ONLY.
 REFER TO CHAPTER 4 FOR DETAILS.
        XRA
                A
                        :CLEAR A.
        OUT
                        ; ACCESS DATA DIRECTION
0003 AF
                       XRA
                                       :CLEAR A.
0004 D301
                       OUT
                               1
                                       :ACCESS DATA DIRECTION
0006 D303
                       OUT
                               3
                                       :REGISTERS.
0008 D305
                       OUT
                               7
000A D307
                       OUT
                                       ;SET AS INPUT.
000C B300
                               0
                       OUT
                               900E,A
000E 3EC0
                       IVM
                               2
0010 B302
                       OUT
0012 3EFF
                       IVM
                               A,377Q
0014 D304
                       DUT
                               4
                                       :SET AS OUTPUT.
0016 D306
                       OUT
0018 3E34
                       IVH
                               A,640
001A D301
                       OUT
                                       :ACCESS DATA REGISTERS.
001C B301
                       OUT
001E D303
                       OUT
                               3
0020 D305
                       OUT
                               5
0022 D307
                       OUT
                IN THE FOLLOWING THE STATUS OF THE
                THE MODE SELETION AUITCH IS DETECTED.
                IF THE NODE IS SELECTED AS TESTING, A
                2KHZ PULSE IS SENT TO EITHER ONE OF THE
                TRANSMITTING TRANSDUCERS. OTHERWISE
                COMPUTATION SEQUENCE IS RESUMED.
                       IN
                              2
                                       ; CHECK TESTING?
0024 DB02
               SNP
0026 E610
                       AHI
                              0208
                              FREE
                                       ; JUMP IF NOT TESTING.
0028 C24000
                       JNZ
002B 3E80
                       IVM
                              A,200Q
                                      COTHERWISE TURN ON PULSE.
                              2
002D D302
                       OUT
                                       :TURN OFF PULSE.
002F 3E00
                       IUK
                              A,0
0031 D302-
                       OUI
                              - 2
0033 210000
                              H,O
                      LXI
0036 7D
              DL
                      HOV
                                       ; DELAY FOR 1/2 MSEC.
                              A,L
0037 23
                       INX
                              Н
0038 FE22
                       CPI
                              420
003A C23600
                       JNZ
                              DL
003B C32400
                       JMP
                              SNP
               START COMPUTATION SEQUENCE.
                ROUTINES ARE SPECIFIED AS FOLLOWING,
```

```
; RST5= RESTORE (H,L> INTO (A)
                ; FACC= FLOATING POINT ACCUMULATOR.
                  FMUL= FLOATING POINT MULTIPLICATION.
                 FSUB= FLOATING POINT SUBTRACTION
                  FDIV= FLOATING POINT DIVISION.
                . FADD= FLOATING POINT ADDITION.
                 FOUT = CONVERT FLOATING POINT NUMBERS INTO BCD.
                 WHERE (A) IS CONTENTS OF A REGISTER.
                ; (H,L) IS LOCATION DEFINED BY REGISTER H AND L.
0040 3E3C
                FREE
                        HVI
                                 A,740
                                         ; CLEAR COUNTER.
0042 D303
                         OUT
                                 3
0044 3E34
                        MVI
                                 A,648
0046 D303
                        DUT
                                 3
0048 3E80
                        MUI
                                 A,2000 ;SEND PULSE TO
004A D302
                        CUT
                                 2
                                         :ULTRASONIC PINGER.
004C AF
                        XRA
                                 A
004B D302
                        OUT
                                 2
004F CD8205
                        CALL
                                 NOTYET
                                          GO CONVERTING ROUTINE.
0052 210A08
                        LXI
                                 H.D1
                                          :STORE 1ST CHORD.
0055 CDCF03
                        CALL
                                 RST
0058 3E3C
                                 A,74Q
                        MVI
                                          CLEAR COUNTER.
005A D303
                                 3
                        OUT
005C 3E34
                        MVI.
                                 A,64Q
005E D303
                        OUT
                                 A,3000
                                        ;SEND 2ND PULSE.
0040 3EC0
                        IVM.
                                 2
0062 D302
                        OUT
                        XRA.
0064 AF
                                 A.
0065 B302
                        CUT
                                 2
0067 CD8205
                        CALL
                                 NOTYET
006A 210E08
                                          :STORE 2ND CHORD.
                        LXI
                                 H,D2
006D CDCF03
                        CALL
                                 RST6
0070 210A08
                                 H,D1
                        LXI
                                          :RESTORE D1.
0073 CD8803
                        CALL
                                 RST5
0076 213108
                        LXI
                                 H, FACC
                                          ;GET DI SQUARE.
0079 CD8001
                        CALL
                                 FMUL
007C 210A08
                        LXI
                                 H,D1
007F CDCF03
                        CALL
                                 RST6
                                          ;STORE DI SQUARE.
0082 210E08
                        LXI
                                 H,D2
0085 CD8803
                        CALL
                                 RST5
0088 213108
                        LXI
                                 H. FACC
008B CD8001
                        CALL
                                 FMUL
                                         GET D2 SQUARE.
008E 210E08
                                 H,D2
                        LXI
0091 CDCF03
                        CALL
                                 RST6
                                         ;STORE D2 SQUARE.
0094 210A08
                        LXI
                                 H.D1
0097 CD9CO4
                        CALL
                                 FADD
                                         :ADD D1^2 AND D2^2.
009A 218004
                        LXI
                                 H.DC
009D CD6005
                        CALL
                                 FSUB
                                         ;SUBTRACT BY D^2.
00A0 211908
                        LXI
                                 H, VA
                        CALL
00A3 CBCF03
                                 RST6
                                         :STORE IT.
00A6 213108
                        LXI
                                 H, FACC
```

; RST6= STORE (A) INTO <H.L>

```
00A9 CB8001
                         CALL
                                  FMUL
                                           ;SQUARE IT.
                         LXI
00AC 218804
                                  H,F4
                                           ;6ET 4.
                         CALL
                                  FDIV
                                           DIVIDED BY 4.
OOAF CDEFO3
00B2 210A08
                         LXI
                                  H,D1
                         CALL
                                  FDIV
                                           :DIVIDED BY D1 2.
00B5 CDEF03
00B8 210E08
                         LXI
                                  H,D2
00BB CDEF03
                         CALL
                                  FDIV
                                           :DIVIDED BY D2^2.
                         LXI
                                  H, VASQ
QOBE 211D08
OOC1 CDCF03
                         CALL
                                  RST6
                                           ;STORE IT.
00C4 218C04
                         LXI
                                  H,FI
00C7 CD8803
                         CALL
                                           ;1 TO FACC
                                  RST5
QQCA 211D08
                         LXI
                                  H. VASQ
                                  FSUB
OOCD CD6005
                         CALL
                                           ;GET 1-(VASQ).
00D0 210E08
                         LXI .
                                  H,D2
00D3 CD8001
                         CALL
                                  FMUL
                                           ;MULTIPLIED BY D2^2.
                         LXI
                                  H,Y2SQ
00D6 212508
                                  RST6 __ ;STORE IT.
QOD9 CDCF03
                         .CALL
                         LXI
                                  H,VA
                                           ;TO RESTORE VA.
00DC 211908
                         CALL
00DF CD8803
                                  RST5
00E2 219004
                         LXI
                                  H,F2
                                           GET 2.
                         CALL
                                 FDIV
00E5 CDEF03
                                           GET VA/2.
                         LXI
                                  H, VAHF
00E8 212108
                         CALL
                                 . RST6
OOEB CDCF03
                                           STORE IT.
                         LXI
                                  H,D2
00EE 210E08
                                           ;RESTORE D2^2.
                         CALL.
00F1 CD8803
                                  RS15
00F4 212108
                         LXI
                                 H,VANF .
                         CALL
                                  FSUB
                                           ;GET D2^2-VAHF.
QOF7 CD6005
                         LXI
00FA 213108
                                  H, FACC
00FD CD8001
                         CALL
                                  FMUL
                                           ;SQUARE IT.
0100 218804
                         LXI
                                  H,F4
                         CALL
                                  FDIV
0103 CDEF03
                                           DIVIDED BY 4.
                         LXI
                                  H.Y2SQ
0106 212508
                         CALL
                                  FDIV
                                           ;DIVIDED BY Y2SQ.
0109 CDEF03
                         LXI
010C 212908
                                  H,YCSQ
                         CALL
                                  RST
010F CDCF03
                                           ;STORE IT.
                         LXI
                                  H,D1
0112 210A08
                         CALL
                                  RST5
                                           ; RESTORE D1^2.
0115 CD8803
                         LXI
0118 218804
                                  H,F4
                                           ;DIVIDED BY 4.
011B CDEF03
                         CALL
                                  FDIV
                         LXI
011E 212908
                                  H,YCSQ
0121 CD9CO4
                         CALL
                                  FADD
                                           ; ADD IT WITH YCSQ.
                         LXI
                                  H.PHI
0124 219404
                         CALL
                                  FMUL
0127 CD8001
                                           ; MULTIPLIED BY PHI.
                         LXI
                                  H,STRIN
012A 214508
                         CALL
                                  FOUT
                                          :CONVERT RESULT INTO BCD.
012D CDF701
0130 CD3601
                         CALL
                                  DISP
                                           DISPLAY RESULT ON LEDS.
0133 C32400
                         JHP
                                  SNP
                                           :GO BACK FOR NEXT.
```

: "<u>.</u> 2 "

```
: THE RESULTS FROM FOUT IS ASSUMED TO BE LESS THAN
                 : 200 CH SQUARE AND GREATER THAN 20 CH SQUARE.
                DISP
                         EQU
0136 =
0136 3E34
                         HVI
                                 A.640
0138 D307
                         OUT
                                 7
                                          :CLEAR 1ST DIGIT.
                                          ; NUMER OF DIGIT INTO C.
                         IVH
                                 C.5
013A 0E05
                                                  ;SET DISPLAY BUFFER.
                                  D,OUTSTR
0130 111208
                         LXI
                                 H,STRIN ;LOCATE BCD RESULT.
 013F 214508
                         LXI
0142 7E
                 REPT
                         YOM
                                 A.M
0143 23
                         INX
                                 Н
                                          ; CHECK IF IS A BCD.
 0144 FE0A
                         CPI
                                  10
                                  REPT
                                          ; IF NOT SKIP IT.
                         JP
 0146 F24201
                                 D
                                          :OTHERWISE INTO BUFFER.
 0149 12
                         STAX
                         INX
                                          :INCREMENT BUFFER.
 014A 13
                                  D
                                          ;LESS ONE DIGIT.
014B OD
                         DCR
                                          GO REPEAT IF MORE DIGIT.
                         JNZ
                                  REPT
014C_C24201
                                                  :LOCATE BUFFER.
                         LXI
                                 H. OUTSTR
 014F-211208
                         VOK
0152 7E
                                  A.H.
 0153 FE02
                         CPI
                                  2
                                          ;CHECK IF >100.0
0155 F25E01
                         JP
                                  SKIP
                                          :SKIP IF NOT.
0158 3E3C
                         MVI
                                  A.740
                                          :DTHERWISE GO LIT
015A B307
                         DUT
                                 7
                                          ;IST DIGIT TO 1.
015C 23
                         INX
                                 Н
015B 7E
                         YOM
                                 A.A
                                          SET NEXT.
015E 23
                SKIP
                         INX
                                 Н
015F 07
                         RLC
                                          ROTATE TO LEFT.
0160 07
                         RLC
0161 07
                         RLC
0162 07
                         RLC
0163 E6F0
                         ANI
                                OFOH-
                                          STRIP TO 4 BITS.
0165 86
                                 ň -.,
                         ADD .
                                          ; MERGE WITH NEXT BCD.
0166 D306
                         OUT
                                 6
                                          :BISPLAY IT(2 DIGITS).
0168 23
                                 Н
                         INX
0169 7E
                         VON
                                          :GET JRD BCD.
                                 A.M
016A 07
                         RLC
                                          ; ROTATE IT.
016B 07
                         RLC
0160 07
                         RLC
016D 07
                        RLC
016E 07
                        RLC
016F E6F0
                        'ANI
                                 OFOH .
0171 23
                         INX
                                 Н
0172 86
                        ADD
                                 H
                                          ; MERGE WITH LAST BCD.
0173 D304
                        OUT
                                 4
                                          :DISPLAY AT LAST 2 DIGITS.
0175 210000
                        LXI
                                 H.O
                                          :DELAY FOR 1 SEC.
0178 7C
                DL1
                        HOV
                                 A,H
0179 23
                        INX
                                 H
017A FE7F
                        CPI
                                 1770
017C C27801
                        JNZ
                                 DL1
017F C9
                        RET
                                         FRETURN.
```

```
FLOATING POINT MULTIPLICATION SUBROUTINE.
                  FLOATING POINT MULTIPLY THE NUMBER AT (H,L) TO
                   THE FACC.
                 EMUL
                        _EDIL
0180 CB6303
                          CALL
                                   FTEST
0183 C8
                          RZ
                          INX
0184 23
                                   H
0185 7E
                          MOV
                                  A,H
                          DCX
                                  Н
0186 2B
0187 B7
                          ORA.
                                  Α
0188 CA8803
                          JZ
                                  RST5
                 FMUL1
018B 113108
                          LXI
                                  D,FACC
018E 1A
                         LDAX
                                  D
018F CD7103
                          CALL
                                  FEXP'
0192 47
                         VOM
                                  B,A
0193 7E
                         VOM
                                   A.M
0194 CD7103
                         CALL
                                  FEXP
0197 80
                         ADD
                                  B
0198 CDD703
                                  FOVUN
                          CALL
                                  127
019B E67F
                         ANI
019D 47
                          MOV
                                  B,A
019E 1A
                         LDAX
                                  D
019F AE
                         XRA
                                  M
01A0 E480
                         ANI
                                  128
01A2 B0
                                  B
                         ORA
01A3 12
                                 - D
                         STAX
01A4 113E08
                         LXI
                                  D,FTEMP+9
01A7 0603
                        · MVI
                                  B,3.
01A9 23
                         INX
                                  Н
01AA CD7803
                                  COPYH
                         CALL
01AB 213508
                         LXI
                                  H,FTEMP
01B0 0606
                         IVM
                                  B,6
01B2 CDE703
                         CALL
                                  ZERON
01B5 113208
                         LXI
                                  B.FACC+1
01B8 0603
                         MUI
                                  B,3
01BA 1A
                FHUL5
                                  D
                         LDAX
01BB 77
                         VOK
                                  M,A
01BC AF
                         XRA
                                  A
01BD 12
                         STAX
                                  D
01BE 13
                                  Ď
                         XKI
01BF 23
                         INX
                                  Н
0100 05
                         DCR
                                  В
01C1 C2BA01
                         JNZ
                                  FMUL5
                         NUI
                                  C.24_
01C4_0E18
                                  B,3
0166 0603
                FNUL6
                         HVI
01C8 213E08
                                  H,FTEMP+9
                         LXI
OICB AF
                         XRA
                                  Α
01CC 7E
                FHULZ
                         KOV
                                  A,H
OICD IF
                         RAR
01CE 77
                         VOM
                                  M,A
```

1. 1

```
01CF 23
                         INX
                                 Н
 01D0 05
                         DCR
                                 B .
 01B1 C2CC01
                         JNZ
                                 FXUL7
 01D4 D2E201
                         JNC
                                 FMUL8
 01D7 113708
                         LXI
                                 D.FTEMP+2
 01DA 213D08
                         LXI
                                 H.FTEMP+8
 01DD 0606
                         IVH
                                 B.6
 01DF CD7D03
                         CALL
                                 FADDT
 01E2 0606
                 FXUL8
                         IVM
                                 B.6
 01E4 213D08
                                 H.FTEMP+8
                         LXI
 01E7 AF
                         XRA.
                                 Α .
 01E8 7E
                 FNUL?
                         VOM
                                 A,N
 01E9 17
                         RAL
 01EA 77
                         HOV
                                 M,A
 01EB 2B
                         DCX
                                 Н
                         DCR
 01EC 05
                                 B
 01ED C2E801
                         JNZ
                                 FMUL9
 01F0 0B
                         DCR
                                 C
 01F1 C2C601
                         JNZ.
                                 FMUL6
 01F4 C39003
                         JNP
                                 FNORM
                 : FLOATING POINT DUTPUT ROUTINE.
                 ; THIS SUBROUTINE CONVERTS A NUMBERS IN THE FLOATING
                 ; POINT ACCUMULATOR TO A BCD STRING SUITABLE FOR DISPLAY.
 01F7 =
                FOUT
                        EBU
                                . $
 01F7 413408
                                D,FACC+3
                        LXI.
 01FA 1A
                        LDAX
                                 B
 01FB F607
                                 7
                        ORI
 01FD 12
                        STAX
                                 D
 01FE CD6303
                        CALL
                                 FTEST
                                M, '
 0201 3620
                        HVI
_ 0203 F20802_
                         JP_
                                 FOUTA ...
                                H, '-'
 0206 362D
                        MVI
 0208 23
                FOUTO
                        INX
                                H 🗀
 0209 C21202
                        JNZ
                                 FOUT2
 020C 3630
                        IVM
                                 M, '0'
 020E 23
                        INX
                                 Н
 020F 3620
                FOUT1
                        IVH
 0211 C9
                        RET
 0212 3A3108
                FOUT2
                        LDA
                                 FACC
 0215 CD7103
                        CALL
                                FEXP
 0218 C21D02
                        JNZ .
                                FOUTV
 021B 3E80
                        HUI '
                                 A,128
 021D E680
                FOUTV
                        ANI
                                 128
021F 320908
                        STA
                                DEXP
0222 E5
                        PUSH
                                Н
0223 3A3108
                FOUT3
                        LDA
                                FACC
0226 CD7103
                        CALL
                                FEXP
0229 FE01
                        CPI
                                1
022B F24402
                        JP
                                FOUT6
```

| 022E 210908 FOUT4 LXI H,D<br>0231 34 INR H   |   |
|--|---|
|  | EXP                                     |
|  |   |
|  | EN                                      |
| 0232 217C04 LXI H,T  |   |
| 0235 F23E02 JP F0U   | TS                                      |
| 0238 CD8001 CALL FHU   | IL                                      |
| 023B C32302 JNP FOU  |   |
|  | -                                       |
| 023E CDEF03 FOUTS CALL FDI   |   |
| 0241 C32302 JNP FOU  | 173                                     |
| 0244 FE05 FOUT& CPI 5  |   |
| 0246 F22E02 JP FOU   | TT A                                    |
|  |   |
|  | TEMP                                    |
| 024C CDCF03 CALL RST   | 6                                       |
| 024F 3A3108 LDA FAC  | C                                       |
| 0252 CD7103 CALL FEX   |   |
|  |   |
| 0255 0E06 NVI C,6  | L                                       |
| 0257 CDAE02 CALL FOU   | TB                                      |
| 025A FE0A CPI 10   |   |
| 025C FA6802 JN FDU   | TIL                                     |
|  |   |
|  | TEMP                                    |
| 0262 CD8803 CALL RST   | <b>'5</b> .                             |
| 0265 C32E02 JMP FOU  | IT4                                     |
| 0268 CD9C02 FOUTU CALL FOU   |   |
|  | 17                                      |
| 026B AF FOUT7 XRA A  |   |
| 026C 323108 STA FAC  | E                                       |
| 026F CD9202 CALL FOU   | TB                                      |
|  | TEMP                                    |
|  |   |
| 0275 CDCF03 CALL: RST  |   |
|  |   |
| 0278 CD9202 CALL FOU   | 118                                     |
|  |   |
| 027B CD9202 CALL FOU   | 18                                      |
| 027B CD9202 CALL FOU<br>027E 113408 LXI D,F  | T8<br>ACC+3                             |
| 027B CD9202         CALL FOU           027E 113408         LXI D,F           0281 213808         LXI H,F   | ACC+3<br>TEMP+3                         |
| 027B CD9202 CALL FOU<br>027E 113408 LXI D,F  | ACC+3<br>TEMP+3                         |
| 027B CD9202         CALL FOU           027E 113408         LXI D,F           0281 213808         LXI H,F   | ACC+3<br>TEMP+3                         |
| 027B CD9202         CALL FOU           027E 113408         LXI D,F           0281 213808         LXI H,F           0284 0604         MVI B,4           0286 CD7D03         CALL FAD  | TEMP+3                                  |
| 027B CD9202         CALL FOU           027E 113408         LXI D,F           0281 213808         LXI H,F           0284 0604         MUI B,4           0286 CD7D03         CALL FAD           0289 CD9C02         CALL FOU   | TEMP+3  BT T9                           |
| 027B CD9202         CALL FOU           027E 113408         LXI D,F           0281 213808         LXI H,F           0284 0604         MVI B,4           0286 CD7D03         CALL FAD           0287 CD9C02         CALL FOU           028C C26B02         JNZ FOU   | ACC+3<br>TEMP+3<br>BT<br>T9             |
| 027B CD9202         CALL FOU           027E 113408         LXI D,F           0281 213808         LXI H,F           0284 0604         MUI B,4           0286 CD7D03         CALL FAD           0289 CD9C02         CALL FOU   | ACC+3<br>TEMP+3<br>BT<br>T9             |
| 027B CD9202         CALL FOU           027E 113408         LXI D,F           0281 213808         LXI H,F           0284 0604         MVI B,4           0286 CD7D03         CALL FAD           0287 CD9C02         CALL FOU           028C C26B02         JNZ FOU           028F C3C502         JMP FOU   | ACC+3<br>TEMP+3<br>BT<br>T9<br>T7<br>TH |
| 027B CD9202         CALL         FOU           027E 113408         LXI         D,F           0281 213808         LXI         H,F           0284 0604         MUI         B,4           0286 CD7D03         CALL         FAD           0287 CD9C02         CALL         FOU           028C C26B02         JNZ         FOU           028F C3C502         JMP         FOU           0292 213408         FOUT8         LXI         H,F   | ACC+3 TEMP+3 DT T9 T7 TH ACC+3          |
| 027B CD9202         CALL FOU           027E 113408         LXI D,F           0281 213808         LXI H,F           0284 0604         MVI B,4           0286 CD7D03         CALL FAD           0287 CD9C02         CALL FOU           028C C26B02         JNZ FOU           028F C3C502         JMP FOU           0292 213408         FOUT8 LXI H,F           0295 54         MOV D,H   | ACC+3 TEMP+3 DT T9 T7 TH ACC+3          |
| 027B CD9202         CALL         FOU           027E 113408         LXI         D,F           0281 213808         LXI         H,F           0284 0604         MVI         B,4           0286 CD7D03         CALL         FAD           0287 CD9C02         CALL         FOU           028C C26B02         JNZ         FOU           028F C3C502         JMP         FOU           0292 213408         FOUT8         LXI         H,Fi           0295 54         MOV         D,H           0296 5D         MOV         E,L  | ACC+3 TEMP+3 DT T9 T7 TH ACC+3          |
| 027B CD9202         CALL FOU           027E 113408         LXI D,F           0281 213808         LXI H,F           0284 0604         MVI B,4           0286 CD7D03         CALL FAD           0287 CD9C02         CALL FOU           028C C26B02         JNZ FOU           028F C3C502         JMP FOU           0292 213408         FOUT8 LXI H,F           0295 54         MOV D,H   | ACC+3 TEMP+3 DT T9 T7 TH ACC+3          |
| 027B CD9202         CALL         FOU           027E 113408         LXI         D,F           0281 213808         LXI         H,F           0284 0604         MVI         B,4           0286 CD7D03         CALL         FAD           0287 CD9C02         CALL         FOU           028C C26B02         JNZ         FOU           028F C3C502         JMP         FOU           0292 213408         FOUTS         LXI         H,F           0295 54         MOV         D,H           0296 5D         MOV         E,L           0297 0604         NVI         B,4   | ACC+3 TEMP+3 BT T9 T7 TH ACC+3          |
| 027B CD9202         CALL         FOU           027E 113408         LXI         D,F           0281 213808         LXI         H,F           0284 0604         MVI         B,4           0286 CD7D03         CALL         FAD           0287 CD9C02         CALL         FOU           028C C26B02         JNZ         FOU           029F C3C502         JMP         FOU           0292 213408         FOUTS         LXI         H,F           0295 54         MOV         D,H           0297 0604         MVI         B,4           0299 C37D03         JMP         FADI  | ACC+3 TEMP+3 BT T9 T7 TH ACC+3          |
| 027B CD9202         CALL         FOU           027E 113408         LXI         D,F           0281 213808         LXI         H,F           0284 0604         MVI         B,4           0286 CD7D03         CALL         FAD           0287 CD9C02         CALL         FOU           028C C26B02         JNZ         FOU           029E C3C502         JMP         FOU           0292 213408         FOUT8         LXI         H,Fi           0295 54         MOV         D,H           0296 5B         MOV         E,L           0297 0604         NVI         B,4           0299 C37D03         JMP         FAD           029C F600         FOUT9         ORI         O  | ACC+3 TEMP+3 BT T9 T7 TH ACC+3          |
| 027B CD9202         CALL         FOU           027E 113408         LXI         D,F           0281 213808         LXI         H,F           0284 0604         MVI         B,4           0286 CD7D03         CALL         FAD           0287 CD9C02         CALL         FOU           028C C26B02         JNZ         FOU           029E C3C502         JMP         FOU           0292 213408         FOUT8         LXI         H,F           0295 54         MOV         D,H           0296 5D         MOV         E,L           0297 0604         NVI         B,4           0299 C37D03         JMP         FADI           029C F600         FOUT9         ORI         O           029E E1         POP         H  | ACC+3 TEMP+3 BT T9 T7 TH ACC+3          |
| 027B CD9202         CALL         FOU           027E 113408         LXI         D,F           0281 213808         LXI         H,F           0284 0604         MVI         B,4           0286 CD7D03         CALL         FAD           0287 CD9C02         CALL         FOU           028C C26B02         JNZ         FOU           029E C3C502         JMP         FOU           0292 213408         FOUT8         LXI         H,Fi           0295 54         MOV         D,H           0296 5B         MOV         E,L           0297 0604         NVI         B,4           0299 C37D03         JMP         FAD           029C F600         FOUT9         ORI         O  | ACC+3 TEMP+3 BT T9 T7 TH ACC+3          |
| 027B CD9202         CALL         FOU           027E 113408         LXI         D,F           0281 213808         LXI         H,F           0284 0604         MVI         B,4           0286 CD7D03         CALL         FAD           0287 CD9C02         CALL         FOU           028C C26B02         JNZ         FOU           029E C3C502         JMP         FOU           0292 213408         FOUT8         LXI         H,F           0295 54         MOV         D,H           0296 5B         MOV         E,L           0297 0604         MVI         B,4           0297 C37D03         JMP         FAD           029C F600         FOUT9         ORI         O           029E E1         POP         H           029F E3         XTHL         VIII   | ACC+3 TEMP+3 BT T9 T7 TH ACC+3          |
| 027B CD9202         CALL         FOU           027E 113408         LXI         D,F           0281 213808         LXI         H,F           0284 0604         MVI         B,4           0286 CD7D03         CALL         FAD           0289 CD9C02         CALL         FOU           028C C26B02         JNZ         FOU           029E C3C502         JMP         FOU           0292 213408         FOUT8         LXI         H,F           0295 54         MOV         D,H           0296 5D         MOV         E,L           0297 0604         MVI         B,4           0299 C37D03         JMP         FAD           029C F600         FOUT9         ORI         O           029E E1         POP         H           029F E3         XTHL         MOV         M,A  | ACC+3 TEMP+3 BT T9 T7 TH ACC+3          |
| 027B CD9202         CALL         FOU           027E 113408         LXI         D,F           0281 213808         LXI         H,F           0284 0604         MVI         B,4           0286 CD7D03         CALL         FAD           0287 CD9C02         CALL         FOU           0286 C26B02         JNZ         FOU           028F C3C502         JMP         FOU           0292 213408         FOUT8         LXI         H,F           0295 54         MOV         D,H           0296 5B         MOV         E,L           0297 0604         HVI         B,4           0299 C37D03         JMP         FADI           029E E1         POP         H           029F E3         XTHL         NA           02A0 77         MSV         H,A           02A1 23         INX         H  | ACC+3 TEMP+3 BT T9 T7 TH ACC+3          |
| 027B CD9202         CALL         FOU           027E 113408         LXI         D,F           0281 213808         LXI         H,F           0284 0604         MVI         B,4           0286 CD7D03         CALL         FAD           0287 CD9C02         CALL         FOU           028C C26B02         JNZ         FOU           029E C3C502         JMP         FOU           0292 213408         FOUT8         LXI         H,F           0295 54         MOV         D,H           0296 5B         MOV         E,L           0297 0604         HVI         B,4           0299 C37D03         JMP         FAD           029E E1         POP         H           029F E3         XTHL         NOV         M,A           02A0 77         MOV         M,A         INX         H           02A1 23         INX         H         MOV         A,C  | ACC+3 TEMP+3 BT T9 T7 TH ACC+3          |
| 027B CD9202         CALL         FOU           027E 113408         LXI         D,F           0281 213808         LXI         H,F           0284 0604         MVI         B,4           0286 CD7D03         CALL         FAD           0287 CD9C02         CALL         FOU           028C C26B02         JNZ         FOU           028F C3C502         JMP         FOU           0292 213408         FOUT8         LXI         H,F           0295 54         MOV         D,H           0296 5D         MOV         E,L           0297 0604         MVI         B,4           0299 C37D03         JMP         FADI           029C F600         FOUT9         ORI         O           029E E1         POP         H           029F E3         XTHL         O           02A0 77         MBV         M,A           02A1 23         INX         H           02A2 79         MOV         A,C           02A3 FE06         CPI         6   | ACC+3 TEMP+3 BT T9 T7 TH ACC+3          |
| 027B CD9202         CALL         FOU           027E 113408         LXI         D,F           0281 213808         LXI         H,F           0284 0604         MVI         B,4           0286 CD7D03         CALL         FAD           0287 CD9C02         CALL         FOU           028C C26B02         JNZ         FOU           028F C3C502         JMP         FOU           0292 213408         FOUT8         LXI         H,F           0295 54         MOV         D,H           0296 5D         MOV         E,L           0297 0604         MVI         B,4           0299 C37D03         JMP         FADI           029C F600         FOUT9         ORI         O           029E E1         POP         H           029F E3         XTHL         O           02A0 77         MBV         M,A           02A1 23         INX         H           02A2 79         MOV         A,C           02A3 FE06         CPI         6   | ACC+3 TEMP+3 BT T9 T7 TH ACC+3          |
| 027B CD9202         CALL         FOU           027E 113408         LXI         D,F           0281 213808         LXI         H,F           0284 0604         MVI         B,4           0286 CD7D03         CALL         FAD           0287 CD9C02         CALL         FOU           028C C26B02         JNZ         FOU           028F C3C502         JMP         FOU           0292 213408         FOUT8         LXI         H,F           0295 54         MOV         D,H           0295 54         MOV         E,L           0297 0604         MVI         B,4           0297 0604         MVI         B,4           0297 F600         FOUT9         ORI         O           029E E1         POP         H           029F E3         XTHL         NOV         A,C           02A0 77         MBV         M,A           02A1 23         INX         H           02A2 79         MOV         A,C           02A3 FE06         CPI         6           02A5 C2AB02         JNZ         FOUT                                 | ACC+3 TEMP+3 BT T9 T7 TH ACC+3          |
| 027B CD9202         CALL         FOU           027E 113408         LXI         D,F           0281 213808         LXI         H,F           0284 0604         MVI         B,4           0286 CD7D03         CALL         FAD           0289 CD9C02         CALL         FOU           028C C26B02         JNZ         FOU           028F C3C502         JMP         FOU           0292 213408         FOUT8         LXI         H,F           0295 54         MOV         D,H           0295 54         MOV         E,L           0297 0604         MVI         B,4           0297 037D03         JMP         FAD           029C F600         FOUT9         ORI         O           029E E1         POP         H           029F E3         XTHL         O           02A0 77         MOV         A,C           02A1 23         INX         H           02A2 79         MOV         A,C           02A3 FE06         CPI         6           02A5 C2AB02         JNZ         FOUT           02A8 362E         HVI         M,* | ACC+3 TEMP+3 BT T9 T7 TH ACC+3          |
| 027B CD9202         CALL         FOU           027E 113408         LXI         D,F           0281 213808         LXI         H,F           0284 0604         MVI         B,4           0286 CD7D03         CALL         FAD           0287 CD9C02         CALL         FOU           028C C26B02         JNZ         FOU           028F C3C502         JMP         FOU           0292 213408         FOUT8         LXI         H,F           0295 54         MOV         D,H           0295 54         MOV         E,L           0297 0604         MVI         B,4           0297 0604         MVI         B,4           0297 F600         FOUT9         ORI         O           029E E1         POP         H           029F E3         XTHL         NOV         A,C           02A0 77         MBV         M,A           02A1 23         INX         H           02A2 79         MOV         A,C           02A3 FE06         CPI         6           02A5 C2AB02         JNZ         FOUT                                 | ACC+3 TEMP+3 BT T9 T7 TH ACC+3          |

| OZAC OD            |        | DCR   | C                |
|--------------------|--------|-------|------------------|
| 02AD E9            | FOUTB  | PCHL  | E A              |
| 02AE 5F<br>02AF AF | FUVID  | MOV   | E,A<br>A         |
|                    |        | XRA   | FACC             |
| 02B0 323108        | EOUTO  | STA   |                  |
| 0283 213408        | FOUTC  | LXI   | H,FACC+3         |
| 02B6 0604          | EQUED. | XVI   | B,4              |
| 02B8 7E            | FOUTD  | MOV   | A,M              |
| 02B9 17            |        | RAL   | 14 A             |
| 02BA 77            | •      | - MOV | M,A              |
| 02BB 2B            |        | DCX   | H                |
| 02BC 05            |        | DCR   | B                |
| 02BD C2B802        | •      | JNZ   | FOUTD            |
| 02C0 1B            |        | DCR   | E                |
| 02C1 C2B302        |        | JNZ   | FOUTC            |
| 02C4 C9            |        | RET   |                  |
| 02C5 E1            | FOUTH  | POP   | H                |
| 02C4 3645          |        | MVI   | H, 'E'           |
| 0208 23            |        | INX   | H                |
| 02C9 3A0908        |        | LDA   | DEXP             |
| 02CC_362B          |        | _NUI_ | \t, <u>_</u> _t_ |
| 02CE 57            |        | HOV   | D,A              |
| 02CF B7            |        | ORA.  | A.               |
| 02D0 F2DC02        |        | JP    | FOUTI            |
| 02D3 362D          |        | HUI   | H, ^-^           |
| 02D5 E67F          |        | ANI   | 127.             |
| 02D7 2F            |        | CHA   |                  |
| 02D8 3C            |        | INR   | A                |
| 02D9 57            |        | VOM   | B,A              |
| 02DA 2F            |        | CMA   |                  |
| 02DB 3C            |        | INR   | A                |
| 02DC 23            | FOUTI  | INX   | <u>H</u>         |
| 02DD E5            |        | PUSH  | H                |
| 02DE 1EFF          |        | MVI   | E,OFFH           |
| 02E0 1C            | FOUTJ  | INR   | E                |
| 02E1 D60A          |        | SUI   | 10               |
| 02E3 F2E002        |        | JP.   | FOUTJ            |
| 02E6 C60A          |        | ADI   | 10               |
| 02E8 47            |        | VOM   | B,A              |
| 02E9 7B            |        | HOU   | A,E              |
| 02EA CD9C02        |        | CALL  | FOUT9            |
| 02EB 78            |        | KOU   | A,B              |
| 02EE CD9C02        |        | CALL  | FOUT9            |
| 02F1 E1            |        | POP   | H                |
| 02F2 3620          |        | IVK   | H, ' '           |
| 02F4 7A            |        | MOV   | A,D              |
| 02F5 B7            |        | ORA   | A                |
| 02F6 F2FF02        | •      | JP    | FOUTK            |
| 02F9 FEFE          | -      | CPI   | OFEH ;-2         |
| 02FB D8            | •      | RC    | , –              |
| 02FC C30203        |        | JMP   | FOUTL            |
|                    |        |       |                  |

| 02FF FE06              | FOUTK  | CPI   | 6                  |     |
|------------------------|--------|-------|--------------------|-----|
| 0301 D0                |        | RNC   |                    |     |
| 0302 4F                | FOUTL  | MOV   | C,A                |     |
| 0303 0605              |        | MVI   | B,5                |     |
| 0305 3620              | FOUTH  | HVI   | M, ′ ′             |     |
| 0307 2B                |        | DCX   | H                  | •   |
| 0308 05                |        | DCR   | B                  |     |
| 0309 C20503            |        | JNZ   | FOUTH              |     |
| 030C EB                |        | XCHG  |                    |     |
| _ 030D ZB              |        | -MOV- | A,E                |     |
| 030E D605              |        | SUI   | 5                  |     |
| 0310 6F                |        | HOU   | L,A                |     |
| 0311 7A                |        | HOV   | A,D                |     |
| 0312 BE00              |        | SBI   | <b>O</b> .         |     |
| 0314 67                |        | HOV . | H,A                |     |
| 0315 79                |        | MOV   | A,C                |     |
| 0316 B7                |        | DRA   | A                  |     |
| 0317 CA2803            |        | JZ    | FOUTO              |     |
| 031A FA3D03            |        | JN:   | FOUTR              |     |
| 031D-46                | FOUTN- | MOV   | B,M                |     |
| 031E 23                |        | INX   | Н                  |     |
| 031F 7E                |        | YOU   | A,M.               | •   |
| 0320 70                |        | YOU   | М,В -              |     |
| 0321 2B                |        | DCX.  | Н                  |     |
| 0322 77                |        | NOV   | H,A                |     |
| 0323 23                |        | INX   | H                  |     |
| 0324 OB                |        | DCR   | C                  |     |
| 0325 C21D03            |        | JNZ   | FOUTN              | · . |
| 0328 EB                | FOUTO  | XCH6  |                    |     |
| 0329 7E                | FOUTP  | MOV   | a,M                |     |
| 032A FE30              |        | CPI   | -0-                |     |
| 032C C33503            |        | JHP   | FOUTO              |     |
| 032F 3620              |        | HVI   | M <sub>2</sub> ( ) |     |
| 0331 2B                |        | DCX   | Н                  |     |
| 0332 C32903            |        | JNP   | FOUTP              |     |
| 0335 FE2E              | FOUTO  | CPI   | <u> </u>           |     |
| 0337 23                |        | INX   | Н                  |     |
| 0338 C0                |        | RNZ   |                    |     |
| 0339 2B                |        | DCX   | H                  |     |
| 033A 3630              |        | HVI   | M, 101             |     |
| 033C C9                |        | RET   |                    |     |
| 033D FEFF              | FOUTR  | CPI   | 255                |     |
| 033F C24B03            |        | JNZ   | FOUTS              |     |
| 0342 2B                | •      | DCX   | H                  |     |
| 0343 7E                |        | HOV   | A,N                |     |
| 0344 362E              |        | MUI   | H, ''              |     |
| 0346 23                | •      | IXX   | H                  |     |
| 0347 77                |        | HOV   | M,A                |     |
| 0348 C32803<br>0348 2B | EDUTO' | JMP   | FOUTO              |     |
| 034C 7E                | FOUTS  | DCX   | <u> </u>           |     |
| 034D 3630              |        | MOV   | A, H               |     |
| VOTU 303V              |        | MUI   | n, '0'             |     |
|                        |        |       |                    |     |

```
034F 23
                     INX
0350 77
                     MOV
                             H.A
0351 62
                     VOM
                             H.D
0352 6B
                     VOK
                             L,E
0353 0606
                     IVM
                             B,6
0355 2B
              FOUTT
                     DCX
                             Н
0356 7E
                     VOK
                             A,X
0357 23
                     INX
                             Н
0358 77
                     YOM
                             N.A
0359 2B
                     DCX
                             Н
035A 05
                     DCR
                             В
035B C25503
                             FOUTT
                     JNZ
035E 362E
                             H, '..'
                     IVM
0360 C32803
                     JNP
                             FOUTO
              ; TEST THE SIGN OF THE NUMBER IN THE FACC.
0363 =
              FTEST
                     EQU
                             $
                             FACC+1
0363 3A3208
                     LDA
0366 B7
                     ORA
0367 C8
                     RZ
0368 3A3108
                             FACC
                     LDA
036B F67F
                     ORI
                             127
036D 3A3108
                             FACC
                     LDA
0370 C9
                     RET
              $P$中国的一个,我们的一个,我们的一个,我们的一个,我们的一个,我们的一个,我们的一个,我们的一个,我们的一个,我们的一个,我们的一个,我们的一个,我们的一个
              : EXPAND EXPONENT INTO 8 BINARY BITS.
0371 =
              FEXP
                     EQU
0371 E67F
                             127
                     ANI
0373 C640
                             64
                     ADI
0375 EE40
                     XRI
                             64
0377 C9
                     RET
              ; HOVE THE STRING FROM (H,L) TO (D,E) COUNT IN B.
0378 =
              COPYH
                     EQU
0378 EB
                     XCHG
0379 CDC603
                             COPYD
                     CALL
                     EXH6
037C C9
                     RET
              ; ADD TWO HULTIPRECISION NUMBERS (D,E) & (H,L).
037B =
              FADDT
                     EQU
037D AF
                     XRA
037E 1A
             FAD1
                     LBAX
                             D
037F 8E
                     ADC
                             M
0380 12
                   · STAX
                             D
0381 1B
                     DCX
                             D
0382 2B
                     DCX
                            H
```

```
0383 05
                      DCR
0384 C27E03
                      JNZ
                             FAD1
0387 C9
                      RET
              ; LOAD THE FLOATING POINT ACCUMULATOR WITH THE 4 BYTES
              ; AT (H,L).
0388 =
              RS15
                      EQU
0388 113108
                      LXI
                             D, FACC
038B 0604
                      IVK
                             B, 4
038D C37803
                      JMP
                             COPYH
              ; NORHALIZE THE FLOATING ACCUMULATOR. THAT IS,
              : THE FIRST BIT MUST BE SIGNIFICANT.
0390 =
              FNORM
                      EQU
0390 213408
                             H.FACC+3
                      LXI
0393 7E
                      MOV
                             A,N
0394 2B
                      DCX
                             Н
0395 B6
                      ORA-
                             N-
0396 2B
                      DCX
0397 B6
                      ORA
                             M
0398 C29E03
                      JNZ.
                             FNRM1
039B 2B
                      DCX
                             Н
039C 77
                      VOM
                             M,A
039D C9
                      RET
039E 7E
              FNRM1
                      MOV
                             A,M
039F B7 -
                      DRA-
                             A٠
03A0 FA6303
                      JM:
                             FTEST
03A3 23
                      INX.
                             H
03A4 23
                      INX
                             H
03A5 0603
                      HUI
                             B,3
                             A
03A7 AF
                      XRA
03A8 7E
              FNRN2
                      YOK
                             A,M
03A9 17
                     RAL
03AA 77
                      VOK
                             M.A.
03AB.2B.
                      DCX
                             H_
03AC 05
                      DER
                             B
03AD C2A803
                             FNRM2
                      JNZ
03B0 7E
                     HOU.
                             A,M
03B1 CD7103
                     CALL
                             FEXP
03B4 3B
                      DCR
0385 FEB0
                     CPI
                             128
0387 CA0000
                     JZ
                             OVERR
038A E67F
                     ANI
                             127
03BC 47
                     YOK
                             B,A
03BD 7E
                     MOV
                             A,N
03BE E480
                     ANI
                             128
03C0 B0
                     ORA
                             В
03C1 77
                     YDK
                             M,A
03C2 23
                     INX
                             Н
```

```
41.
```

```
03C3 C39E03
                         FNRN1
                   JMP
            : HOVE THE STRING FROM (D.E) TO (H.L) COUNT IN B.
0306 =
            COPYD
                  EQU
                         $ .
03C6 1A
                   LDAX
                         B
0307 77
                  YOH
                         M,A
03C8 23
                   INX
                         H
0309 13
                   IXX
                         D
03CA 05
                  DCR
03CB C2C603
                   JNZ
                         COPYD
O3CE C9
                  RET
0000 =
            OVERR
                   EQU
            ; STORE THE FLOATING POINT ACCUMULATOR AT (H,L)
03CF =
            RST6
                  EQU
                         $
                  LXI
03CF 113108
                         D.FACC
03D2 0604
                  MUI
                         B.4
03D4 C3C6O3
                   JMP
                         COPYD
            ; TEST EXPONENT FOR OVERFLOW OR UNDERFLOW.
            FOVUN-
03D7 =
                  EQU
                         $
03D7 B7
                  ORA
                         A
                   JP.
03D8_F2F103
                         FRU1
03DB FEC1
                  CPI.
                         193
03DD DO
                  RNC
03DE C30000
                  JHP
                        OVERR
03E1 FE40
            FOV1
                  CPI
                         64
03E3 D8
                  RC
03E4 C30000
                  JMP
                         OVERR
            ; HOVES A STRING OF BINARY HZEROS. COUNT IN B.
03E7 =
            ZERON
                  EQU
                         $
03E7 3600
                  HVI
                         M.O
03E9 23
                  INX
                         Н
03EA 05
                  DCR
                         B
03EB C2E703
                  JNZ
                         ZEROM
03EE C9
                  RET
            : FLOATING POINT DIVISION SUBROUTINE.
            ; FLOATING POINT DIVIDE THE NUMBER AT (H.L) INTO ..
            ;THE FACE.
03EF =
            FDIV
                  EOU
                         $
03EF CD6303
                  CALL
                         FIEST
03F2 C8
                  RZ
03F3 23
                  INX
                         Н
                         A,H
03F4 7E
                  MOV
```

| 03F5 2B              |         | DCX    | Н         |
|----------------------|---------|--------|-----------|
| 03F6 B7              |         | ORA    | Á         |
| 03F7 CA0000          |         | JZ     | OVERR     |
| 03FA -7E             |         | VON    | A,M       |
| 03FB CD7103          |         | CALL   | FEXP      |
| 03FE 47              |         | MOV    | B,A       |
| 03FF 113108          |         |        | •         |
|                      |         | LXI    | D, FACC   |
| 0402 1A              |         | LDAX   | . D .     |
| 0403 CD7103          |         | CALL   | FEXP      |
| 0406 90              |         | SUB    | В         |
| 0407 3C              |         | INR    | Α         |
| 0408 CDD703          |         | CALL   | FOVUN     |
| 040B E67F            |         | ANI    | 127       |
| 040D 47              |         | VOV    | B,A       |
| 040E 1A              |         | LDAX   | D         |
| 040F AE              |         | XRA    | М         |
| 0410 E680            | •       | ANI    | 128       |
| 0412 B0              |         | ORA    | В         |
| 0413 12              | •       | STAX   | D .       |
| 0414 E5              |         | PUSH   | H         |
| 0415 13              |         | INX    | · D       |
| 0416 213508          |         | LXI    | H,FTEMP   |
| 0419 3600            |         | HVI    | н,о       |
| 0417 3300<br>041B 23 |         | INX    | H         |
| 041C 0603            |         |        |           |
|                      | FD707   | MVI    | B,3       |
| 041E 1A              | FDIV3   | LDAX   | D         |
| 041F 77              | •       | HOV    | M,A       |
| 0420 AF              |         | XRA    | A         |
| 0421 12              |         | STAX   | D         |
| 0422 23              | •       | INX    | H         |
| 0423 13              | . •     | INX    | D         |
| 0424 05              |         | , DCP, | В         |
| 0425 C21E04          |         | JNZ    | FDIV3     |
| 0428 D1              |         | POP    | D         |
| 0429 0603            | *       | IVH    | B,3       |
| 042B 13              |         | INX    | D         |
| 042C 3600            |         | HVI    | N,0       |
| 042E 23              | •       | INX    | H         |
| 042F CDC603          |         | CALL   | COPYD     |
| 0432 0E18            |         | MVI    | C,24      |
| 0434 113808          | FDIV5   | LXI    | D,FTEMP+3 |
| 0437 213008          |         | LXI    | H,FTEMP+7 |
| 043A 0604            |         | HVI    | B,4       |
| 043C CD7104          |         | CALL   | FSUBT     |
| 043F D24E04          |         | JNC    | FDIV6     |
| 0442 113808          |         | LXI    | D,FTEMP+3 |
| 0445 213008          |         | LXI    | H,FTENP+7 |
| •                    |         |        | •         |
|                      |         | MUI    | B,4       |
| 044A CD7D03          | · · ·   | CALL   | FADDT     |
| 044D 37              | PR-711. | STC    |           |
| 044E 3F              | FDIV6   | CMC    |           |

```
044F 0603
                        MUI
                                B,3
                                H,FACC+3
0451 213408
                        LXI
               EDIU7
                        אמע
0454 7E
                                A_H_
                        RAL
0455 17
0456 77
                        YOM
                                M,A
                        DCX
                                H
0457 28
                        DCR
0458 05
                                В
0459 C25404
                        JNZ
                                FDIV7
                        XRA
045C AF
                                A
045B 0604
                        MVI
                                B.4 .
045F 21380B
                        LXI
                                H.FTEMP+3
0462 7E
               FDIV8
                        VON
                                A.A
0463 17
                        RAL
0464 77
                        MOV
                                H,A
0465 28
                        DCX
                                Н
                        DCR
                                B
0466 05
0467 C26204
                        JNZ
                                FDIV8
046A OD
                        DCR
046B C23404
                        JNZ
                                FDIV5
046E C39003
                        JMP
                                FNORM
                ; SUBTRACT THE TWO HULTIPRESION NUMBERS (D.E) 1 (H,L)
0471 =
                FSUBT
                        EQU
                        XRA
0471 AF
                                A
0472 1A
               FSB1
                        LDAX
                                D
                        SBB
                                M
0473 9E
                        STAX
0474 12
                                B
0475 1B
                        DCX
                                D
                        DCX
0476 28
                                Н
0477 05
                        DCR
                                B
0478 C27204
                        JNZ
                                FSB1
047B C9
                        RET
                : FOLLOWING ARE CONSTANTS TO BE USED IN THE COMPUTATION
                 SEQUENCES. ALL THE CONSTANTS ARE ARRANGED IN FLOATING
                ; POINT FORMAT.
047C 04
               TEN
                        DB
                                04H
                                         ; INTEGER 10.
                        DB
                                OAOH
047D A0
                        DB
047E 00
                                0
047F 00
                        DB
                                0
               DC
                        DB
                                        ;SQUARE OF 2.5(CM).
0480 03
                                03H
                        DB
                                0C8H
0481 C8
                        DB
                                OOH
0482 00
                        DB
                                OOH
0483 00
               TUD
0484 02
                        DB
                                02H
                                         ;INTEGER 2.
                        DB
                                80H
0485 80
                        DB
0486 00
                        DB
0487 00
                                0
               F4
                        DB
                                03H
0488 03
                                         ; INTEGER 4.
0489 80
                        DB
                                80H
```

```
048A 00
                        DB
                                0
048B 00
                                0
                        DB
                F1
048C 01
                        DB
                                1
                                         :INTEGER 1.
048D 80
                        DB
                                80H
048E 00
                        DB
                                00
048F 00
                                00
                        DB
0490 02
                F2
                        DB
                                2
0491 80
                        DB
                                80H
0492 00
                        DB
                                0
0493 00
                        DB
                                0
0494 02
                PHI
                                2
                        DB
                                         :PHI=3.1415926
0495 C9
                                0C9H
                        DB
0496 OF
                        DB
                                OFH
0497 B7
                        DB
                                0D7H
                SPEED
0498 7A
                                         :SPEED OF ULTRASOUND.
                        DB
                                7AH
0499 BC
                        DB
                                8CH
                                         ;=0.15CH/USEC.
049A A4
                        DB
                                0A4H
049B 9C
                        DB
                                9CH
                : FLOATING POINT ADDITION SUBROUTINE.
                ; FLOATING POINT ADD THE NUMBER AT (H,L) TO FACC
049C. =
                FADD
                        EQU
                              . $
049C 23
                                H
                        INX
049D 7E
                        YOM
                               A,H
049E B7
                        ORA
                                A
                                FTEST
049F CA6303
                        JZ -
04A2 2B
                        DCX
                                H
04A3 CD63O3
                        CALL
                                FTEST
04A6 CA8803
                        JZ
                                RST5
04A9 113108
                        LXI
                                D,FACC
04AC 1A
                        LDAX
                                D
04AB CB7103
                        CALL
                                FEXP
04B0 47
                        VOK
                                B,A
04B1 7E
                        VOK
                                M.A
0482 CD7103
                        CALL
                                FEXP
04B5 4F
                        VOK
                                C,A
04B6 90
                        SUB
04B7 CACCO4
                                FADD4
                        JZ
                        JP.
                                FADD3
Q4BA F2BEQ4
04BD 2F
                        CMA
04BE 3C
                        INR
                                A
04BF FE18
               FADD3
                        CPI
                                24
04C1 DACC04
                        JC
                                FADD4
04C4 78
                        YOM
                                A,B
04C5 91
                        SUB
                                C
04C6 F26303
                                FTEST
                        J۶
04C9 C38803
                        JMP
                                RST5
04CC F5
               FADD4
                       PUSH
                                PSW
04CD C5
                       PUSH
                                В
04CE 113508
                        LXI
                                D.FTEMP
```

|   | 04D1 | 0604       |        | MVI  | B,4        |
|---|------|------------|--------|------|------------|
|   | 04D3 | CD7803     |        | CALL | COPYH      |
|   | 0406 | C1         |        | POP  | В          |
| , | 04D7 | F.1        |        | POP  | PSW        |
|   |      | CAFF04     |        | JZ   | FADD9      |
|   | ,    | 213608     |        | LXI  | H,FTEMP+1  |
|   | 04DE |            |        | PUSH | PSU        |
|   | 04DF |            |        | YON  | A, B       |
|   | 04E0 |            | •      | SUB  | c ·        |
|   | -    | F2F104     |        | JP   | FADD6      |
|   |      | 213108     |        | LXI  | H.FACC     |
|   | 04E7 |            |        | HOU  | A,C        |
|   |      | E67F       |        | ANI  | 127        |
|   | 04EA |            |        | MOV  | C,A        |
|   | 04EB |            |        | HOU  | A,M        |
|   |      | E680       |        | ANI  | 128        |
|   | 04EE |            |        | ORA  | C          |
|   | 04EF | 77         |        | VOM  | N,A        |
|   | 04F0 |            | •      | INX  | H.         |
|   | 04F1 | FI         | FADD6  | POP  | PSV        |
|   | 04F2 | 4F         |        | HOV  | C,A        |
|   | 04F3 | 0603       | FADD7  | MVI  | B,3        |
|   | 04F5 | AF         |        | XRA- | A          |
|   | 04F6 | <b>E</b> 5 | **     | PUSH | H          |
|   | 04F7 | CD7905     |        | CALL | FSHFT      |
|   | 04FA | E1         |        | POP  | H          |
|   | 04FB | OD         |        | DCR  | C          |
|   | 04FC | C2F304     |        | JNZ. | FADD7      |
|   | 04FF | ÷          | _FADD2 | EOU  | <b>_\$</b> |
|   | 04FF | 113108     |        | LXI  | B,FACC     |
|   | 0502 | 213508     |        | LXI  | H,FTEMP    |
|   | 0505 | 1 A        |        | LDAX | D          |
|   | 0506 | AE         |        | XRA  | М.         |
|   | 0507 | FA3205     |        | JM   | FADDA      |
|   | 050A | 113408     | •      | LXI  | D,FACC+3   |
|   |      | 213808     | •      | FXI  | H,FTEMP+3  |
|   |      | 0603       | •      | NVI  | B,3        |
|   |      | CB7B03     |        | CALL | FADBT.     |
|   |      | D26303     |        |      | FTEST      |
|   |      | 213108     |        | LXI  | H,FACC     |
|   | 0518 |            |        | MOV  | A,M        |
|   |      | E680       |        | ANI  | 128        |
|   | 051E |            |        | MOV. | B,A        |
|   | 051F |            | •      | HOV  | A, H       |
|   |      | CD7103     |        | CALL | FEXP       |
|   | 0523 |            |        | INR  | A          |
|   |      | E67F       |        | ANI  | 127        |
|   | 0526 |            |        | ORA  | В          |
|   | 0527 |            | .•     | MOV  | N,A        |
|   | 0528 | 23         |        | INX  | Н          |
|   |      |            |        |      |            |

```
0529 37
                       STC
052A 0603
                       IVM
                               B,3
052C CB7905
                               FSHFT
                       CALL
052F C36303
                       JNP
                               FTEST
0532 =
               FADDA
                       EQU
0532 213908
                       LXI
                               H.FTEMP+4
0535 0604
                       IVM
                               B.4
0537 CDC603
                       CALL
                               COPYD
               FADDB
                               D,FACC+3
053A 113408
                       LXI
053D 213808
                       LXI
                               H,FTEMP+3
                               B,3
0540 0603
                       HUI
0542 CD7104
                       CALL
                               FSUBT
0545 D29003
                       JNC
                               FNORM
0548 113108
                       LXI
                               D.FACC
054B 21350B
                       LXI
                               H,FTEMP
054E 0608
                       IVM
                               B,8
0550 CD7803
                       CALL
                               COPYH
0553 113108
                       LXI
                               D.FACC
0556-213508
                       LXI
                              H. ETENP
                       YON.
0559 7E
                               A,N
055A EE80
                       XRI
                               128
055C. 12
                       STAX
                               D
055D C33A05
                       JMP
                               FADDB
                ******************************
               : FLOATING POINT SUBTRACTION SUBROUTINE.
                FLOATING POINT SUBTRACT THE NUMBER AT (H,L)
               FROM THE FACC.
0560 =
                       EQU
               FSUB
0560 23
                       INX
                               H
0561 7E
                      MOV
                              A,M
0562 B7
                       DRA
                              Α
0543 CA4303
                       JZ
                              FTEST
0566 2B
                      DCX
0567 113508
                      LXI
                              D.FTEMP
056A 0605
                      IVM
                              B,5
056C CD7803
                      CALL
                              COPYH
056F 213508
                      LXI
                              H.FTEMP
0572 7E
                      MOV
                               A,M
0573 EE80
                      XRI
                              128
0575 77
                      VON
                              N,A
0576 C39C04
                       JNP
                              FADD
                INCREMENTING SHIFT RIGHT.
0579 =
               FSHFT
                      EQU
0579 7E
                      HOV
                              A,M
057A 1F
                      RAR
057B 77
                      VOM.
                              M.A
```

```
057C 23
                       INX
                               H
057D 05
                       DCR
                               Ř
057E C27905
                       JNZ
                               FSHFT
0581 C9
                       RET
               SUBROUTINE "NOTYET".
                 THIS ROUTINE TAKES A BINARY NUMBER
                 FROM THE COUNTER (THE COUNTER WHICH
                 HAS THE TRAVERSE TIME INFORMATION)
                 AND CONVERT IT INTO CORESPONDING
                CHORD LENGTH.
                 THE RESULT IS IN FLOATING POINT
               : FORMAT AND IS NORWALIZED.
0582 =
               NOTYET
                      EQU
                               $
0582 DB02
                       IN
                               2
0584 E620
                       ANI
                               0400
                                       :CHECK IF COUNTER READY?
0586 C28205.
                       JNZ
                              NOTYET
                                       :GO BACK IF NOT.
0589 DB00
                       IN
                                       COTHERWISE TAKE 1ST BYTE.
                       YOM
058B 5F
                               E.A
                                       :STORE IT.
058C 07
                       RLC
                                       :ROTATE IT.
058D 07
                     RLC
058E 07
                       RLC
058F 07
                       RLC
                      ANI
                               OFOH.
                                       ;STRIP INTO 4 BITS.
0590 E&FO
                                       STORE INTO 3RD BYTE OF FACC.
0592 323308
                       STA
                               FACC+2
0595 7B
                      VOK
                               A,E
                                       :RESTORE IST BYTE.
                                       ; ROTATE IT.
0596 OF
                      RRC
0597 OF
                       RRC
0598 OF
                       RRC
0599 OF
                       RRC
059A E60F
                       ANI
                               OFH
                                       STRIP TO 4 BITS.
                                       :STORE IN E.
059C 5F
                               E,A
                       MOV
                               2
059B BB02
                       IN
                                       :GET NEXT HIGHER BYTE.
                       RLC
059F 07
                                       ; ROTATE.
05A0 07
                       RLC
05A1 07
                       RLC
05A2 07
                       RLC
                               OFOH
05A3 E6F0
                       ANI
                                       :STRIP TO 4.
05A5 83
                       ADD
                                       : MERGE.
05A6 323208
                       STA
                               FACC+1
                                       :STORE INTO FACC.
                             A.OCH
                                       ; HOVE '12' INTO EXPONENT
05A9 3E0C
                      MVI.
                               FACC
05AB 323108
                       STA
                                       :OF FACC.
                       IVA
                                       CLEAR LAST BYTE OF
05AE 3E00
                               A,0
05B0 323408
                       STA
                               FACC+3
                                       ;FACC.
05B3 CD9003
                       CALL
                               FNORM
                                       ; NORMALIZED IT.
05B6 219804
                       LXI
                               H.SPEED ; MULTIPLIED COUNTS BY SPEED
05B9 CD8001
                      CALL
                               FHUL
                                       :TO GET CHORD LENGTH.
05BC C9
                       RET
                                       RETURN.
```

|      | F<br>Skrik & ikiskiski | kale datelesk date | koleskok & skolesk k | E :le:koled | ::::::::::::::::::::::::::::::::::::: | 2.注生体治体体水学体际体等体体 |
|------|------------------------|--------------------|----------------------|-------------|---------------------------------------|------------------|
|      | <b>y</b>               | JING ARE           |                      |             |                                       |                  |
|      | ;                      |                    |                      |             | ,                                     |                  |
| 0800 |                        | ORG                | 4000Q                |             |                                       |                  |
| 0800 | TENP1                  | DS                 | .4                   |             |                                       | •                |
| 0804 | TEMP2                  | DS                 | 4                    |             |                                       |                  |
| 0808 | TEM                    | DS                 | 1                    |             | •                                     | •                |
| 0809 | DEXP                   | DS                 | 1                    |             |                                       |                  |
| 080A | D1                     | DS                 | 4                    |             |                                       | •                |
| 080E | 12                     | DS                 | 4                    |             |                                       |                  |
| 0812 | OUTSTR                 | DS                 | ż                    |             |                                       |                  |
| 0819 | VA                     | DS                 | 4                    |             | •                                     |                  |
|      |                        |                    | 7                    |             |                                       |                  |
| 081D | VASQ                   | DS                 | 4                    |             |                                       |                  |
| 0821 | VAHF                   | DS                 | 4                    |             |                                       |                  |
| 0825 | Y2\$0                  | DS                 | 4                    |             |                                       |                  |
| 0829 | YCSQ                   | DS                 | 4                    |             | ··                                    |                  |
| 082D | VB                     | DS                 | 4                    | •           | ,                                     |                  |
| 0831 | FACC                   | DS                 | 4                    |             |                                       |                  |
| 0835 | FTEHP                  | DS                 | 12                   |             |                                       |                  |
| 0841 | MULPR                  | DS                 | <u> </u>             |             |                                       |                  |
| -    |                        |                    | 7                    |             |                                       |                  |
| 0845 | STRIN                  | DS                 | 20                   |             |                                       |                  |
| 0859 |                        | END                |                      |             |                                       |                  |

#### REFERENCES:

- 1. Paulev, P. E., F. Neumann, S. L. Nielson and N. Keiding. "Strain gage versus water plethysmography; Description of simplified system and analysis of differences and accuracy." Med. Biol. Eng. 12, pp 437-445, 1974.
- 2. Whitney, R. J. "The measurement of volume changes in human limbs." J. Physiol. 12, pp. 1-27, 1953.
- 3. Sigdell, J. E., "Venous occlusion plethysmography." Biomedical Engineering, August and September 1973.
- 4. Lafferty, J. F., "A feasibility study of limb volume measuring system. Final report NAS9-12749, September 1974 and Final report NAS9-14392, October 1977.
- 5. Berry, C. A. and A. D. Catterson. "Pre-Gemini medical prediction versus Gemini flight results." Gemini Summary Conference, NASA SP-138, 1967.
- 6. Dietlein, L. F. and W. V. Judy. "Cardiovascular Conditioning." A Review of Medical Results of Gemini 7 and Related Flights, NASA TM X-60589, pp. 1-35, 1966.
- 7. Hoffler, G. W. and R. L. Johnson. "Apollo flight crew cardiovascular evaluations." In <u>Biomedical Results of Apollo</u>, NASA, 1975, pp. 227-264.
- 8. Hoffler, G. W., R. L. Johnson and R. A. Wolthuis. "Apollo space crew cardiovascular evaluations." <u>Aerospace Medicine</u>, 45: pp. 807-820, 1974.
- 9. Johnson, R. L., G. W. Hoffler, A. E. Nicogossian, S. A. Bergman, Jr. and M. M. Jackson. "Lower body negative pressure: Third Manual Skylab Mission." NASA, SP-377, pp. 284-312, 1977.
- 10. Thornton, W. E. and G. W. Hoffler. "Hemodynamic studies of the legs under weightlessness." NASA, SP-377, pp. 324-329, 1977.
- 11. Barnett, R. D., R. J. Gowen, F. L. Zaebst and D. R. Carrol. "A Capacitance Plethysmograph." Seventh Annual Rocky Mountain Bio-Engineering Symposium, 1970.
- 12. Bhagat, P., M. Kadaba, R. Ware and W. Cockerill. "Frequency dependence of acoustic parameters of freshly excised tissues of sprague dawley rats." <u>Ultrasonics</u>, pp. 179-182, July, 1977.
- 13. Kadaba, M. P. "Tissue characterization using ultrasonic transmission and scattering parameters." Ph.D. Dissertation, University of Kentucky, 1978.
- 14. Rushmer, R. F., D. L. Franklin and E. M. Ellis. "Left ventricular dimensions recorded by sonocardiometry." <u>Circulation Research</u>, Volume 4, pp. 684-688, 1956.

- 15. Stegall, H. F., M. B. Kardon, H. L. Stone and V. S. Bishop. "A portable, simple sonomicrometer." J. Appl. Physiol., 23, 289-293, 1967.
- 16. Horwitz, L. D., V. S. Bishop, H. L. Stone and H. F. Stegall. "Continuous measurement of internal left ventricular diameter." J. Appl. Physiol., 24, 738-740.
- 17. Stegall, H. F., J. M. Ofstad, V. E. Simmons and R. F. Rushmer. "An integrated facility for peripheral vascular studies in man." Volume 6, p. 75, Proceeding 17th ACEMB, Cleveland, 1964.
- 18. Bhagat, P. K., J. F. Lafferty and M. P. Kadaba. "An ultrasonic plethysmograph." Aerospace Medical Association Meeting, Washington, D. C., Reprint, 73-74, May 1979.
- 19. Kadaba, M. P., P. K. Bhagat and V. Wu. "Attenuation and backscattering of ultrasound in freshly excised animal tissues." IEEE Trans. BioMed Eng., BME-27(2), 76-83, 1980.
- 20. Bhagat, P. K., J. F. Lafferty, D. Bowman and M. P. Kadaba. "An ultrasonic plethysmograph for space flight applications." Aviat. Space Environ. Med. 51(2), 185-188, 1980.
- 21. Wu, V. C. "A microprocessor based ultrasonic lower limb volume measurement system." M. S. Thesis, University of Kentucky, 1980.
- 22. Levitan, B. M. "Evaluation of plethysmographic systems proposed for use on the space shuttle." Special Report NAS9-14880, 1980 and Aviat. Space Environ. Med. 53(1), 6-10, 1983.
- 23. Bhagat, P. K., M. P. Kadaba, V. N. Gupta and V. Wu. "Microprocessor based system for ultrasonic tissue characterization." Med. Instrum. 14(4), 220-224, 1980.
- 24. Dekker, A. J. "Solid state physics", Prentice-Hall, Sixth Edition, 184-209, 1962.
- 25. Mason, W. P., "Piezoelectric crystals and their application to ultrasonics." D. Van Nostrand Co., pp. 234-288, 1950.
- 26. Hueter, T. F. and Bolt, <u>Ultrasonics</u>, 59-85 and 401-415, John Wiley, 1955.
- 27. Van Trigt, P., B. J. Bauer, C. O. Olsen, J. S. Rankin and A. S. Wechsler. "An improved transducer for measurement of cardiac dimensions with sonomicrometry." Am. J. Physiol., 240 (Heart Circ. Physiol. 9), H664-H669, 1981.
- 28. DeSilets, C. S., J. D. Fraser and G. S. Kino. "The design of efficient broadband piezoelectric transducers." <u>IEEE Trans. Sonics Ultrason.</u>, SU25, 115-125, 1978.